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
EN

GENERAL

Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at service@vanguard.de or contact our field staff before using the product on patients.

General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

- Instructions
- Lists

Device information

This product may only be operated if its safe use is guaranteed. Observe the intended purpose and safety instructions!

Vanguard assumes no liability for damages that arise from improper operation of the Vanguard Irrigated Ablation Catheter NAV ST BW and/or disregard of these instructions for use.


Vanguard AG and the manufacturer of the original product, Biosense Webster, Inc. (33 Technology Drive, Irvine, California 92618, USA), are not affiliated companies; the device is remanufactured solely by Vanguard AG, who markets the device without collaboration by Biosense Webster, Inc. "THERMOCOOL SMARTTOUCH" and "CARTO" are protected brands of the manufacturer of the original product, Biosense Webster Inc. (33 Technology Drive Irvine, California 92618, USA) or one of its affiliated companies.

"Redel" is a registered, protected brand of INTERLEMO HOLDING S.A. (INTERLEMO HOLDING AG / INTERLEMO HOLDING LTD); Chemin des Champs Courbes 28, 1024 Ecublens, Switzerland (CH).

PRODUCT DESCRIPTION

The Vanguard Irrigated Ablation Catheter NAV ST BW is intended for electrophysiological mapping and RF ablation in the heart. The catheter consists of a 7.5 F shaft (2) with a distal 8 F catheter tip (3) and a handle (1). The catheter tip consists of a head electrode (8) with irrigation holes (10) and five ring electrodes (9). Ablation follows with the head electrode. This contains an integrated temperature sensor for temperature monitoring. The head electrode also features an integrated force sensor, which measures the contact force between the head electrode and the heart wall in real time. The head electrodes of the catheters are also equipped with a magnetic sensor. This serves to visualise the position of these catheters with a compatible 3D navigation system. The electrical connection (4) for the RF generator and the 3D navigation system is located proximally on the handle. The application area is irrigated and cooled with heparinised saline solution through the shaft lumen via the luer connection (5) proximal to the handle. This requires the connection of a compatible irrigation pump. Potential recording and stimulation can be performed with the three ring electrodes and the head electrode, which are positioned in the working area.

The distal catheter tip is controlled using the control elements on the handle: curvature control (6) and friction control (7, bidirectional only). The unidirectional catheters BWV and BWX can be deflected in one direction, the bidirectional catheters BWW and BWY in two. Correct positioning of the catheter tip is supported by the torsion-resistant shaft. The types of curvature are listed in the following table. The curvature type is marked on the label.

INTENDED PURPOSE

Indication

The Vanguard Irrigated Ablation Catheter NAV ST BW is indicated in conjunction with compatible accessories for catheter-based electrophysiological mapping (stimulation and recording) and in conjunction with a compatible RF generator for cardiac ablations.

In combination with a compatible CARTO® 3 navigation system, the Vanguard Irrigated Ablation Catheter NAV ST BW provides real-time information on the contact force between the catheter tip and the heart wall as well as on the position of the catheter tip.

Guidelines for use

The medical device may only be used by a trained and experienced specialist physician. All applied medical procedures and techniques are the responsibility of the attending physician. The information contained in these instructions of use must be observed.

Contraindications

The Vanguard Irrigated Ablation Catheter NAV ST BW is contraindicated:

- in patients having undergone ventriculotomy or atriotomy within the last 8 weeks
- in patients with intracardiac thrombus or myxoma
- in patients with prosthetic heart valves
- in coronary vessels
- in patients with active systemic infection
- in the event of transseptal access in patients with interatrial deflections or corrections
- in patients intolerant to anticoagulation therapy (heparin or equivalent alternative)
- when using a sheath or introducer <8.5 Fr

Side effects

A number of severe side effects have been documented for catheter ablation procedures. These include pulmonary embolism, myocardial infarction, strokes, pericardial tamponade and death.

The following complications have also been reported in previous studies or are listed in literature.

In connection with catheterisation / catheterisation procedures:

Vascular bleeding/local haematoma, thrombosis, arteriovenous (AV) fistulae, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, pericardial effusion/tamponade, thrombi, air embolism, arrhythmia and valve injuries, pneumothorax and haemothorax, water retention and lung oedema, hypoxia, pleural effusion, Acute Respiratory Distress Syndrome (ARDS), congestive cardiac insufficiency, aspiration pneumonia, pneumonia, asthma attacks, hypotension, malfunction of the implantable Cardioverter Defibrillator (ICD), anaemia, thrombocytopaenia, disseminated intravascular coagulation, haemorrhagic areas, systemic infections, urinary tract infections, apnoea caused by sedation, CO2 retention caused by sedation with lethargy and cholecystitis.

In connection with RF:

Chest pain/complaints, ventricular tachyarrhythmia, transitory ischaemic attack (TIA), stroke (CVA), complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery section, thromboembolism of the heart, pericarditis, perforation of the heart/pericardial tamponade, valve damage and elevated phosphokinase levels.

Independent of devices or procedures:

Urine retention, transient numbness of the limbs, Parkinson's disease and gastrointestinal haemorrhages. Previously unknown side effects may occur at any time. The distributing company must be notified of any side effects that are not specified. In the event of known side effects, please also refer to appropriate expert literature.

TECHNICAL DATA

Diameter:	7.5F (2.5 mm)
Max. diameter:	8F (2.67 mm)
Usable length:	115 cm ± 1,15 cm
Irrigation port:	Luer ISO 594-2 female Variant A
Sensors in catheter tip:	Temperature sensor (Thermocouple) Magnetic sensor, force sensor
Irrigation holes:	Non SF: lateral on head electrode SF: lateral and distal at the head electrode
Length of head electrode:	3,5 mm
Number of ring electrodes:	5
Distance between electrodes:	(1 – 6 – 2) mm
Rated voltage:	220 V _p

BWV	uni-directional	curvature	D	F	J	–	–
		REF No.	35117	35118	35519	–	–
	non SF	Corresponding REF no. of the OEM	D133601	D133602	D133603	–	–

BWW	bi-directional	curvature	D-D	F-F	J-J	F-J	D-F
		REF No.	35123	35124	35125	35126	35127
	non SF	Corresponding REF no. of the OEM	D132701	D132702	D132703	D132704	D132705

BWX	uni-directional	curvature	D	F	J	–	–
		REF No.	35120	35121	35122	–	–
	SF	Corresponding REF no. of the OEM	D134701	D134702	D134703	–	–

BWY	bi-directional	curvature	D-D	F-F	J-J	F-J	D-F
		REF No.	35128	35129	35130	35131	35132
	SF	Corresponding REF no. of the OEM	D134801	D134802	D134803	D134804	D134805

This product bears a CE mark in accordance with the requirements in the EC directive 93/42/EEC.

All Vanguard Irrigated Ablation Catheter NAV ST BW are as compatible with the accessories (e.g. generator, cable) as the corresponding products of the original manufacturer Biosense Webster (see table). When the catheter is connected to an RF generator that provides an RF output voltage > 220 V_p, the maximum RF output voltage must be set to 220 V_p. Only approved medical devices may be used during the procedure with strict adherence to and understanding of their instructions for use.

SAFETY INSTRUCTIONS

1. Do not apply excessive force if you encounter resistance when advancing or retracting the catheter.
2. Excessive bending or buckling can damage the catheter. Manual bending of the distal end may result in damage to the control mechanism and injury to the patient.
3. Contact force measurement serves information purposes only. The standard handling precautions must be observed.
4. As indicated in the handling instructions, the catheter must be preheated prior to use to achieve a constant, stable condition. Otherwise, there is a risk of zero offset, which could result in incorrect force measurement.
5. The contact force must always be set to zero after insertion into the patient and when changing into the other ventricle. Ensure that the catheter is not in contact with the heart tissue during zero-point adjustment.
6. Contact force measurement can be impaired if the contact force sensor (between the 1st and 2nd ring electrode) comes close to ferrous material, e.g. the shaft of another catheter. In the event of extreme force fluctuations, check the zero point and whether the contact force sensor is located near another catheter shaft. It may be necessary to retract the catheter for inspection.
7. To ensure the function of the contact force sensor, the distal three ring electrodes of the catheter tip must protrude from the introducer sheath.
8. If high lateral forces are exerted during mapping and RF application, the user should observe the contact force and vector display on the CARTO® 3 system to ensure that the contact force remains within the desired range.
9. The temperature values transmitted by the sensor in this catheter are not the actual temperatures of the tissue. The temperature values only represent the values of the irrigated electrode. The temperature sensor serves to monitor the irrigation flow rate. Before RF ablation is activated, a drop in the electrode temperature confirms the beginning of the irrigation of the head electrode with saline solution. The maintenance of the irrigation flow rate during RF ablation is ensured by monitoring the electrode temperature.
10. Do not rely on the electrode temperature to rise during RF ablation to check whether the tissue is heated.
11. The output adjustment procedure described in the instructions for use must be followed. If the output is increased too quickly during ablation, sudden steam development may cause perforation.
12. To avoid overhydration, the patient's fluid balance should be monitored throughout the application. Proceed according to your hospital protocol. In some patients, factors may occur that affect the ability to cope with overhydration. These patients may experience pulmonary oedema or heart failure during or after application. A particular risk exists in patients with congestive cardiac insufficiency, renal insufficiency and in older patients. The risk of overhydration of the patient must be determined before starting the application.
13. There are no studies on the risks of discontinuation of anticoagulation therapy after catheter ablation of atrial fibrillation; for such patients, anticoagulation therapy must be continued in accordance with the "ACC/AHA/ESC Guidelines for the Management of Patients in Atrial Fibrillation".
14. To avoid thromboembolism, intravenous heparin should be used when RF ablation is performed in the left ventricle. Follow the clinical guidelines and general practice for optimal anticoagulation therapy after the procedure.
15. There are no studies on the safety and efficacy of RF ablation for the treatment of atrial fibrillation in patients with significant left ventricular dysfunction, advanced heart failure, significant left ventricular atrial enlargement or structural heart disease.
16. Pacemakers and implantable cardioverters/defibrillators can be impaired or damaged by RF signals. Discharges from intracardiac defibrillators (ICD) during application may result in injury to the patient. The following precautions should therefore be taken during use: a) Provide external stimulation and defibrillation during ablation; b) Deactivate intracardiac defibrillators (ICD) as they may be discharged and injure the patient or damage the patient during ablation; c) Take special care during ablation near permanent stimulation or defibrillation electrodes; and d) Carefully inspect implantable aggregates in all patients following ablation.
17. Following septal accessory pathway ablation, patients are at an increased risk of a complete AV block. A complete AV block requires a permanent pacemaker. A complete AV block may occur after RF ablation.
18. Adequate fluoroscopic visualisation is required for access through the aorta to avoid the risk of catheter placement in the coronary vessels. Ablation in and near arterial coronary vessels has been associated with myocardial infarction and death.
19. Due to the duration and intensity of X-rays for fluoroscopic imaging during catheter ablation procedures, there is a risk of significant exposure to radiation with acute radiation damage and an increased risk of somatic and genetic effects for the patient and laboratory staff. A catheter ablation procedure may only be performed after thorough consideration of the potential radiation exposure and after the implementation of adequate radiation protection measures. The long-term risks of prolonged exposure to X-rays were not determined. For pregnant women and pre-pubescent children, a catheter ablation procedure must therefore be carefully considered.
20. The catheter must not be brought into contact with organic solvents such as alcohol.
21. Do not immerse the handle and cable plug in liquid. This can impair the electrical function.
22. The electrode on the catheter tip must not be rubbed, twisted or pulled during cleaning. This could loosen the electrode from the shaft or damage the force sensor. A significant change in the contact force in the relaxed state after cleaning can indicate damage to the force sensor.
23. Before inserting the catheter into the vascular system, be sure to remove all air from the tubes and catheter as trapped air can cause injury or death.
24. Cardiac ablation procedures may only be performed by physicians who have mastered the techniques of RF frequency catheter ablation. These procedures must be performed in a fully equipped electrophysiological laboratory.
25. No information is available on the long-term risks of RF-induced lesions. The use of the catheter in prepubescent children must be carefully considered.
26. Before inserting and withdrawing the catheter, always straighten the curvature using the control elements.
27. Always ensure continuous irrigation with normal heparinised saline solution to prevent coagulation in and around the head electrode.
28. RF ablation is interrupted if the limit value of the temperature or the impedance is exceeded. In this case, the catheter must be removed from the patient and any clots removed from the distal tip. If the catheter is defective, it must be replaced by a new catheter. Before reinserting the catheter, ensure that the irrigation ports are unblocked.
29. A low output power, a high impedance value or improper functioning of the accessory with normally set parameters can indicate an incorrect application of the neutral electrode(s) or failure of an electrical line. In this case, the current should only be increased after obvious defects, incorrect use of the neutral electrode(s) or the failure of electrical lines have been excluded.
30. It has been proven that ablation systems with irrigation produce larger lesions than RF ablation catheters without irrigation. Special care should be taken when ablating near structures that are sensitive to electrical charges and near thin-walled or other arterial structures.
31. This product is intended for single use only. It may not be reprocessed or reused. Reuse may result in product failure, patient injury and/or transmission of infectious diseases from one patient to another patient or third party.
32. The materials of the catheter are not compatible with magnetic resonance imaging (MRI).
33. The catheter must be advanced under fluoroscopic imaging control. This minimises the risk of heart damage, perforation or tamponade. The compatible catheter positioning and navigation system can be used in conjunction with fluoroscopic imaging.
34. The Vanguard Irrigated Ablation Catheter NAV ST BW emits electrical energy during RF ablation. Improper handling of the catheter or neutral electrode(s) may result in injury to the patient or user. During RF ablation, the patient must not come into contact with earthed metal objects.
35. Electrosurgical applications pose the risk of ignition of flammable gases or other materials. Flammable materials should be kept away from the electrosurgical site.
36. RF ablation generates electromagnetic interference (EMI). This may negatively influence the function of other devices.
37. Electrodes and probes from monitoring and stimulation devices can form paths for high-frequency currents. They should therefore be placed as far away from the ablation site and neutral electrodes as possible to reduce the risk of burns. Protective impedances can reduce the risk of burns. This allows continuous monitoring of the electrocardiogram during the release of RF energy.
38. Prior to RF ablation, ensure that the temperature display of the generator is in working order. The catheter's temperature measurement accuracy is mainly determined by the accuracy of the temperature measurement function of the RF generator used. Observe the respective instructions for use.
39. When performing RF ablation near the phrenic nerve, precautions must be taken to avoid injury to the phrenic nerve. Necessary measures include stimulation to determine the proximity of the catheter tip to the nerve and reduction of RF energy.

HANDLING INSTRUCTIONS

Storage

The device must be stored in dry and light-protected conditions at 5 °C to 25 °C. Information about storage conditions is given on the label on the sterile packaging.

Removal from packaging

- Check that the packaging is not damaged and check the validity of the expiry date. Do not use the device if the packaging is damaged or opened and the expiry date has expired.
- Remove the catheter from its packaging under aseptic conditions and transfer it to a sterile working environment. Do not throw the instrument.
- Check the catheter for damage such as kinks in the shaft, loose electrodes, damaged plugs or defects such as clogged irrigation ports, irrigation openings and lumen, as well as difficulty in adjusting the distal end and insufficient fixation of the distal end. Do not use damaged or inadequate catheters.

Application

- Under consideration of the aseptic rules, access a large central blood vessel through which you can insert the catheter.
- Use the supplied introducer sleeve to advance or retract the catheter through the haemostatic valve of the sheath to prevent damage to the catheter tip. After inserting the catheter, push the introducer sleeve back to the handle.
- To check the compatibility of the catheter with the sheath, pass the catheter through the sheath before inserting it into the patient. Sheaths < 8.5 F are contraindicated.
- Use a compatible cable to connect the catheter to compatible accessories: the CARTO® 3 system, the recording and mapping system and the RF generator. Connect a neutral electrode to the generator to close the circuit. Connect the irrigation pump to the catheter's luer connection. If required, you can use a 3-way stopcock. Connect a bag of normal heparinised saline solution (1 IU heparin/ml) at room temperature to the irrigation pump in accordance with a safe standard procedure.
- Irrigate the catheter and irrigation lines using a safe standard procedure to remove all air bubbles from the catheter and tube and check that all irrigation openings are unblocked.
- The catheter must be continuously irrigated at a continuous flow rate of 2 ml/min.
- Insert the catheter through the puncture site using the supplied introducer sleeve and an introducer sheath of compatible size. Set the parameters on the RF generator depending on the selected catheter type non SF / SF.
- Check the correct positioning of the catheter fluoroscopically.
- For an optimally accurate and stable contact force measurement, the catheter must be warmed up for 2 minutes after connection to the CARTO®-3 system.
- Perform a zero offset of the contact force after inserting the catheter into the patient. The head electrode and all four distal ring electrodes must be outside the sheath to ensure that the force sensor is located in the body. Ensure that the catheter tip does not touch tissue when performing zero offset by checking the position via EMG amplitude with the CARTO® 3 system and the catheter movement under fluoroscopy. Changes in the contact force at the same frequency as the cardiac and respiratory circulation may indicate contact between the tip of the catheter and tissue. If all markers mentioned confirm that there is no contact between catheter tip and tissue, the contact force can be set to zero.
- Repeat the zero offset of the contact force when you move the catheter from one ventricle to another or reinsert it after retraction.
- The distal catheter tip is deflected using the control elements on the handle. With unidirectional catheters, the curvature is set by feeding and retracting the curvature control. With bidirectional catheters, the curvature is set by rotating the curvature control. The degree of the catheter tip's curvature is determined by the strength of the deflection of the control element. The unidirectional catheter is in the neutral position when the control element is retracted. The bidirectional catheter is in the neutral position when the control element is in the middle position.
- Ensure that the RF generator does not display a temperature above 37°C after inserting the catheter into the patient and before delivering RF energy. Before ablation, ensure that the impedance of the circuit is within the typical range.
- Recommended settings for RF ablation and the flow rate:

	Atrial ablation		Ventricular ablation	
Output range	15 W to 30 W		31 W to 50 W	
Application period	30 to 120 seconds		60 to 120 seconds	
Type of irrigation hole	non SF	SF	non SF	SF
Temperature monitoring	< 50 °C	< 40 °C	< 50 °C	< 40 °C
Irrigation flow rate during RF application	17 ml/min	8 ml/min	30 ml/min	15 ml/min

- It is recommended to increase irrigation to the high flow rate 5 seconds before RF ablation and to maintain it for 5 seconds after RF ablation.
- RF energy must not be activated until the increase in the flow rate is confirmed by a temperature drop of at least 2 °C at the electrode tip.
- The temperature at the catheter tip must be constantly monitored during use to ensure adequate irrigation. At a temperature above 40 °C for SF catheters and 50 °C for non SF catheters, RF ablation should be interrupted. The irrigation system must be checked for blockages before renewed RF ablation.
- Starting from the recommended initial output of 15 W to 20 W, the output can be increased in 5 W to 10 W increments after about 15 seconds until a transmural lesion is achieved. Transmural lesions are characterised by an > 80% reduction in unipolar atrial electrogram amplitude or by the occurrence of double potential at equal or lower amplitudes. An output of max. 50 W should be selected if the catheter tip is parallel to the tissue. If the tip of the catheter is at right angles to the tissue, select a maximum output of 35 W. Individual RF ablations should not exceed 120 seconds, during which the catheter can be pulled to the next ablation site. RF ablation can be reapplied at the same site.
- RF ablation can be performed several times at the same site. If RF ablation is interrupted by the generator after exceeding the impedance or temperature limit values, the catheter must be withdrawn and the electrode tip checked for coagulate before it is used for further RF ablations. To remove coagulate, carefully clean the tip with a gauze swab moistened with saline solution. Ensure that all irrigation openings are clear before reinserting the catheter:

Clogged irrigation holes must be cleaned as described below:

- Fill a 1 ml or 2 ml syringe with sterile saline solution and connect it to the catheter (e.g. at the side arm or stopcock).
- Carefully inject the saline solution into the catheter. You should be able to observe an even flow of fluid at the tip of the catheter.
- Repeat steps a) and b) until all irrigation holes are unblocked.
- Irrigate the catheter and irrigation lines using a safe standard procedure to remove all air bubbles from the catheter and tube and check that all irrigation openings are unblocked.
- Reinsert the catheter into the patient.
- After repositioning, carry out a zero offset for the contact force.

NOTE: Outputs > 30 W can be used if no transmural lesions could be produced at lower output levels.

NOTE: When treating isthmus-dependent flutters, output levels above 30 W (up to max. 50 W) should only be selected if no conduction block could be achieved at lower output levels.

NOTE: The temperature displayed on the RF generator only indicates the temperature in the electrode. This is not the temperature of the tissue or that of the interface between the tissue/electrode.

NOTE: In the bidirectional catheter, the degree of curvature fixation can be adjusted all the way to a full lock. This is set with the friction control element. The degree of fixation increases the further the knob is turned clockwise (in plus direction). Delivery follows in non-fixed setting.

Disposal

The product must be disposed of in accordance with the applicable national laws and directives.

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©Vanguard AG 1000001984_REV_C_26-08-2021