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Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at service@vanguard.de or contact our field staff before using the product on patients.

General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

- ✓ Requirements
- Instructions
- Lists
-  Safety instructions

The original version of these instructions for use is in German.

Device Information

This product may only be operated if its safe use is guaranteed. Observe the safety and handling instructions!

VANGUARD assumes no liability for damages that arise from improper operation of the Vanguard Ablation Catheter Cable MTC | OSB and/or disregard of these instructions for use.

Vanguard AG and the manufacturers of the original products, Medtronic, Inc. (710 Medtronic Parkway, Minneapolis, MN 55432, USA) and Osypka AG (Earl-H.-Wood-Straße 1, 79618 Rheinfelden, Germany) are not affiliated companies; the device is remanufactured solely by Vanguard AG, who markets the device without collaboration by Medtronic, Inc. and/or Osypka AG.

“Redel” is a protected brand of the manufacturer INTERLEMO HOLDINGS S.A., Ecublens, CH or of a manufacturer associated with this company.

PRODUCT DESCRIPTION

The Vanguard Ablation Catheter Cable MTC | OSB are flexible, sterile connection cables intended for single use. Plug connectors, which can be designed in several standardized variants, are located at both ends of the product.

INTENDED PURPOSE

Indication

The Vanguard Ablation Catheter Cable MTC/OSB are intended to connect a compatible RF catheter to a compatible RF generator within the scope of an ablation procedure.

Guidelines for use

It is the responsibility of the physician to use a suitable surgical procedure and technology. The procedure described in the instructions for use is for information purposes only. Each physician must apply, supplement or adapt the information in the instructions for use according to his/her medical training and clinical experience. The medical device may only be used by trained personnel in a fully equipped electrophysiology laboratory.

Contraindications

There are no known contraindications for the Vanguard Ablation Catheter Cable MTC | OSB.

TECHNICAL DATA

REF.	Corresponding REF. No. (of the original manufacturer)	Length in [m]	Plug connector Catheter connection	Plug connector System connection	Compatible type of temperature sensor	Output to be set at the RF generator
35590	4825 (Medtronic, Inc.)	1.83	Redel 10-pole, green	Redel 14-pole, black	Thermo-couple (TC)	≤ 50 W
35591	82521 (Ospyka AG)	2.0	Redel 9-pole, blue	Redel 10-pole, red	Thermistor (TR)	≤ 100 W
35592	82514 (Ospyka AG)	1.5	Redel 9-pole, blue	Redel 10-pole, red	Thermistor (TR)	≤ 100 W
35593	81997 (Ospyka AG)	1.5	Redel 10-pole, blue	Lemo 10-pole, metal	Thermo-couple (TC)	≤ 100 W

Other data → see label.

All Vanguard Ablation Catheter Cable MTC | OSB are as compatible with RF generators and RF ablation catheters as the corresponding products of the original manufacturers (see table).

This product bears a CE mark in accordance with the requirements in the EC directive 93/42/EEC and all applicable changes.

SAFETY INSTRUCTIONS

The safety instructions consist of three aspects: Cause, Consequence, Countermeasure. A safety instruction indicates a potential risk to the health or life of persons.

Plug and/or cable comes into contact with liquid!

Impairment with the function of the medical device

- The plug and the cable may not come into contact with liquids.
- The product may not come into contact with organic solvents.

Do not touch the plug contacts with your fingers or thin objects!

Impairment with the function of the medical device

- Ensure that the plug contacts do not come into contact with thin objects. Do not touch the plug contacts.

Lack of caution when handling the Vanguard Ablation Catheter Cable!

Injuries to the patient, the user and/or the OR personnel

- Take great care when handling the Vanguard Ablation Catheter Cable.
- Observe the instructions for use of the Vanguard Ablation Catheter Cable.

Ignition of flammable gases or other substances!

Fire or explosion

- Ensure that flammable and easily inflammable materials are removed from the working area.

Damaged product or product is damaged during use!

Extension of the operating time, impaired RF generator performance

- Carry out a careful visual and functional inspection of the product.
- Never use damaged products.

Sterile packaging damaged or opened!

Infection of the patient due to non-sterile products

- Check the packaging before opening it. The cable may not be used if its sterile packaging is open or damaged.

Damage due to unsuitable storage conditions, use after use-by-date!

Impaired electrical function of the product, incorrect ablation site, death of the patient

- Take note of the relevant labelling on the packaging.

→ Do not use the product once its use-by-date has expired.

Reprocessing, resterilization, reuse and/or technical modification of the cable by the user or operator!

Incorrect energy transfer, uncontrolled release of RF energy, contamination of the patient, infection, death of the patient

→ Vanguard Ablation Catheter Cable MTC | OSB are single use products and intended for single use on one patient, only.

→ Vanguard Ablation Catheter Cable MTC | OSB may not be reprocessed, resterilized, reused or technically modified by the user or operator.

HANDLING INSTRUCTIONS

Handling the cables requires great care to avoid injuries to the patient.

Storage

→ Storage conditions are provided on the label on the product packaging.

Safety measures before use

- The instructions for use of the Vanguard Ablation Catheter Cable MTC | OSB must have been read and understood before use.
- When connecting external devices and using them for specific applications, observe the corresponding instructions for use. Only apply and use authorized medical devices under consideration of the respective instructions for use.
- Any static electricity must be kept away from the vicinity of the cable and its associated configuration. Take particular care to ensure that the operating table and the electronic devices used (e.g. X-rays) are sufficiently and centrally earthed. The cable may not touch the patient and other lines. Electrostatic voltages between patient and user must be equalized before the cable is used.
- A prerequisite for the use of the Vanguard Ablation Catheter Cable MTC | OSB is the appropriate training of medical personnel and the existence of a fully equipped electrophysiology laboratory.

Removal from packaging

- Check the sterilization packaging for integrity and then check the sterilization indicator. The product is sterile if the EO indicator is set to yellow. Products with damaged packaging or non-yellow EO indicator must be disposed of and a new product must be used for the procedure.
- Take the cable out under aseptic conditions and transfer it to a sterile working environment.
- Ensure that the packaging material is disposed of safely.
- Carry out a careful visual and functional inspection of the product. Check for defects such as kinks or damage to the plug. The product may not be used in the event of damage or malfunction.

Application

- ✓ The principles of asepsis must be observed.
- ✓ An emergency power supply must be available to ensure continuous treatment.
- ✓ Defibrillation, emergency pacemaker therapy and acute cardioversion must be possible during the entire treatment.

- ✓ An external pacemaker, defibrillator, intubation kit, oxygen connection and emergency medication must be ready and available throughout the procedure.
- ✓ The procedure must always be performed under X-ray control and permanent ECG monitoring.
- ✓ In order to prevent delays in the process, a sterile substitute should be available during treatment.
- Insert the plugs of the cable into the corresponding connection on the RF generator and into the RF ablation catheter, observing the marking. Check that the plugs are seated correctly.
- Upon completion of the ablation procedure: To disconnect the plug connectors, release the connector lock and pull the plug connector from the respective connection.

NOTE: The Vanguard Ablation Catheter Cable MTC | OSB must be used under sterile conditions.

NOTE: If the product is damaged during use, e.g. kinks in the cable or damage to the pins in the plug, it must be replaced with a new one.

Disposal

The product may be contaminated after use and must therefore be considered as contaminated hospital waste. The cable must be disposed of in accordance with the applicable national laws and directives.

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©Vanguard AG 1000002196 Revision_B_10.03.2021