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GENERAL

Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at service@vanguard.de or contact our field staff before using the product on patients.

General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

- ✓ Requirements
- Instructions
- Lists
-  Safety instructions

The original version of these instructions for use is in German.

Device information

This product may only be operated if its safe use is guaranteed. Observe the notes and safety precautions!

VANGUARD assumes no liability for damages that arise from improper operation of the Vanguard Ablation Catheter uni LT G BTA/ LT BTB/ G BTC/ BTD and/or disregard of these instructions for use.

Vanguard AG and the manufacturer of the original product, VascoMed GmbH (Hertzallee 1, 79589 Binzen, Germany), are not affiliated companies; the device is remanufactured solely by Vanguard AG, who markets the device without the cooperation of VascoMed GmbH and BIOTRONIK SE & Co. KG.

PRODUCT DESCRIPTION

The Vanguard Ablation Catheters uni LT G BTA/ LT BTB/ G BTC/ BTD are catheters with four electrodes and flexible tips. These catheters are intended for electrophysiological mapping of the heart and conduct an alternating current into the tissue via the distal catheter tip for the purpose of ablation. Ablation requires the appropriate equipment (e.g. RF generator, cable, neutral electrode). The catheters are equipped with a torsion-resistant shaft and flexible tip, which contains three ring electrodes and a head electrode. The existing ring electrodes are for discharge and stimulation purposes and the head electrode for ablation purposes. The Vanguard Ablation Catheters uni LT G BTA/ G BTC have gold electrodes and the Vanguard Ablation Catheters uni LT BTB/ BTD have platinum/iridium electrodes. The catheter tip contains a temperature sensor, with which the temperature of the head electrode can be monitored during ablation. The curvature of the tip is controlled by a thumb switch in the hand piece at the proximal end of the catheter. To bend the tip, the thumb switch is pushed in distal direction. When the thumb switch is in the proximal position, the tip is straight. Catheters are available with different curvatures. The designation of the curvature types corresponds to the designations of the original manufacturer VascoMed GmbH: "FullCircle Red", "FullCircle Blue", "FullCircle Green", "FullCircle Black" and "FullCircle Cyan". The level at which the tip is positioned can be rotated by turning the shaft. This

©Vanguard AG 1000001845_REV_D_26-08-2021

allows the practitioner to achieve a specific position of the catheter tip in conjunction with tip flexion and rotation. The catheter is used in combination with suitable recorders and RF generators. For connection, use appropriate patient cables with a suitable plug connection.

Intended purpose

Indication

The Vanguard Ablation Catheters uni LT G BTA/ LT BTB/ G BTC/ BTB are intended for use in conjunction with compatible connecting cables, introducer sheaths, RF generators and suitable devices for intracardiac RF ablation, as well as for recording intracardiac signals and for diagnostic stimulation. The Vanguard Ablation Catheters uni LT G BTA/ LT BTB/ G BTC/ BTB are intended for the treatment of supraventricular, as well as ventricular and atrioventricular tachyarrhythmias or AV reentry tachycardias.

Guidelines for use

It is the responsibility of the physician to use a suitable surgical procedure and technology. The procedure described in the instructions for use is for information purposes only. Each physician must apply, supplement or adapt the information in the instructions for use according to his/her medical training and clinical experience. The medical device may only be used by a trained and experienced cardiologist specialised in electrophysiology.

All Vanguard Ablation Catheters uni LT G BTA/ LT BTB/ G BTC/ BTB are as compatible with the accessories (e.g. generator, cable) as the corresponding products of the original manufacturer VascoMed GmbH.

Contraindications

Use of the Vanguard Ablation Catheter uni LT G BTA/ LT BTB/ G BTC/ BTB is contraindicated in patients with intracardiac mural thrombus; who have undergone ventriculotomy or atriotomy in the previous four weeks; with prosthetic heart valves; with left atrial myxomas or myxomas in the left atrium; with intracardiac patch or baffle (for transseptal access); with a prosthetic aortic valve (for retrograde transaortic access); who have a history of sensitivity to intracardiac patch or baffle (for retrograde transaortic access). intracardiac patch or baffle (for transseptal access); a prosthetic aortic valve (for retrograde transaortic access); a history of sensitivity or allergy to foreign objects; massive cardiac or vascular anomalies (transposition) and vascular obstruction; with unstable angina pectoris; with bacteraemia or sepsis; with acute decompensated congestive heart failure not caused by arrhythmia; with haemorrhagic diathesis; with acute thrombosis in the lower limbs (for femoral catheter access); with atrial thrombi or myxomas in the chambers or target chambers to be passed; acute systemic or local infection; previously diagnosed coronary artery spasm; haemodynamic instability; hypercoagulability; confirmed or suspected obstruction in the vein or artery selected for insertion; physiological or anatomical abnormalities that may cause postoperative complications (e.g., reduced resistance to infection); with severe concomitant diseases resulting in prolonged treatment with significantly increased risk; with damaged tissue (e.g. tissue damaged by radiation); with recurrent metastatic cancer; with active implants (e.g. pacemakers, defibrillators); with bradycardic events; with total AV block and presence of contraindications to transseptal puncture.

Side effects

Irrespective of the proper use of the Vanguard Ablation Catheters uni LT G BTA/ LT BTB/ G BTC/ BTB, the following complications cannot be ruled out during or after treatment: Acute coronary syndrome/occlusion; air embolism; reaction to anaesthesia; stroke, TIA, cerebral lesions; arteriovenous fistulae; atrial fistulae; arterial spasm; cardiac arrest; chest pain and discomfort; cardiac arrhythmias (ventricular and supraventricular tachycardia, atrial fibrillation/flutter, ventricular fibrillation, bradycardia, sinus or AV block/dysfunction,

asystole, extrasystole); diaphragmatic paralysis; thrombosis; embolism; oesophageal fistula; oesophageal ulceration; gastroparesis; haemothorax; haematoma; hypotension/hypertension; hypovolaemia; injury and perforation of nerves, arteries and veins; arterial and venous dissection; myocardial infarction; valvular injury, valvular damage; cardiac injury or perforation/ tamponade; pericarditis, endocarditis, perimyocarditis; phrenic nerve damage; pleural or pericardial effusion; pneumothorax; skin burns; stenosis; thoracic paraesthesia; vasovagal reactions; unretrieved device fragments (UDF); infection; pulmonary embolism; pericardial tamponade; arterial spasm; reversible conduction disturbances due to mechanical contact with the ablation catheter and conduction system; irreversible conduction disturbance which may alter heartbeat; death.

As a consequence of the significant radiation exposure (X-ray exposure), radiation damage can occur. The patient and user are subject to the risk of resulting genetic and somatic damage.

During the procedure, there are risks of ventricular fibrillation caused by spurious currents, as well as arrhythmias caused by pacing in the vulnerable phase or by catheter manipulation. Patients with active implants (e.g. pacemakers) may experience damage to or dislocation of the electrodes as well as disturbances in pacemaker function.

When using the catheter, the following undesirable side effects in addition to those already listed may occur due to catheter breakage, catheter dislocation, contact problems between the catheter and the recording device and an increase in the stimulus threshold: Elevation of phosphokinase levels; subcutaneous haematoma; skeletal muscle and nerve stimulation; lesion of the brachial plexus; lesion of the thoracic duct; lesion of the intracardiac conduction system; pneumothorax during puncture of the subclavian vein; pseudoaneurysms; cardiac thromboembolism; coronary artery spasm; coronary artery thrombosis; coronary artery transection; loss of conduction efficacy (e.g. exit block); cardiac arrest; adverse subclavian artery puncture; venous inflow tract injury.

Adverse reactions are mentioned in literature in connection with cardiac catheterisation and catheter ablation. These are mentioned below in addition to those already listed: Injury to the endocardium up to perforation of the heart wall; vascular haemorrhage, coronary sinus perforation; stenosis of the coronary artery or pulmonary vein; pneumothorax during puncture of the subclavian vein; haematothorax during puncture of the subclavian vein; formation of atriooesophageal fistulas; decompensation as a result of cardiac/renal insufficiency.

Previously unknown side effects may occur at any time. The distributing company must be notified of any side effects that are not specified.

Technical data

Vanguard Ablation Catheter uni LT G BTA/ LT BTB/ G BTC/ BTD

Shaft diameter: 7F
 Usable length: 110 cm
 Electrode distance: 2-5-2 mm
 No. of electrodes: 4
 Temperature sensor: Thermoelement Type T
 Connection: 9 PIN (Redel PRB female)

Vanguard Ablation Catheter uni LT G BTA

REF	Corresponding REF no. of the original	Head	Curvature type	Electrode material
-----	---------------------------------------	------	----------------	--------------------

	product	electrode length		
34977	364539	8 mm	FullCircle Red	Gold
34976	364527		FullCircle Blue	
34978	364535		FullCircle Green	
34971	364524		FullCircle Black	
34979	364531		FullCircle Cyan	

Other data → see label.

Vanguard Ablation Catheter uni LT BTB

REF	Corresponding REF no. of the original product	Head electrode length	Curvature type	Electrode material
34980	364526	8 mm	FullCircle Blue	Platinum/Iridium
34972	364523		FullCircle Black	

Other data → see label.

Vanguard Ablation Catheter uni G BTC

REF	Corresponding REF no. of the original product	Head electrode length	Curvature type	Electrode material
34981	364537	4 mm	FullCircle Red	Gold
34982	364525		FullCircle Blue	
34983	364533		FullCircle Green	
34969	364522		FullCircle Black	
34984	364529		FullCircle Cyan	

Other data → see label.

Vanguard Ablation Catheter uni BTD

REF	Corresponding REF no. of the original product	Head electrode length	Curvature type	Electrode material
34985	364536	4 mm	FullCircle Red	Platinum/Iridium
34986	359038		FullCircle Blue	
34987	364532		FullCircle Green	
34970	359037		FullCircle Black	

Other data → see label.

This product bears a CE mark in accordance with the requirements in the EC directive 93/42/EEC and all applicable changes.

SAFETY INSTRUCTIONS

The safety instructions consist of three aspects: Cause, consequence and counter-measure. A safety instruction indicates a potential risk to the health or life of persons.

-  Proximal plug connection or handle or cable comes into contact with liquid!
Interference with the function of the catheter
→ Prevent the plug, the cable and the handle from coming into contact with liquids.
-  Plug connector damaged!
Interference with electrical safety and function
→ Carry out a careful visual inspection and functional inspection of the plug.
-  Temperature or impedance set point exceeded!
Interruption or influence of the ablation current
→ Bring the catheter tip to a straight position with the aid of the thumb switch and then pull it out and clean coagulation deposits from the tip. Leave the electrode tip in a straight position when cleaning. Avoid rotating or rubbing to prevent damage to the connection point and the electrode tip or loosening of the electrodes.
-  More than 30 ablations conducted each with one minute of energy supply
Hazard due to loss of function or loose electrodes
→ Do not perform more than 30 ablations with an energy supply of over one minute or less each.
-  No anticoagulant is administered!
Application for thromboembolism
→ Administer intravenous heparin into the left ventricle during access for ablation purposes. There is still no consensus concerning the need to use anticoagulants after ablation.
-  The medical personnel have not been adequately trained!
Fault when making intracardiac potential recordings and/or stimulation and ablation.
→ It is a pre-requirement for intracardiac potential recording and/or stimulation that specialist medical personnel are trained sufficiently and that a fully equipped electrophysiology laboratory is present. The procedure itself may only be carried out by a trained and experienced cardiologist specialised in electrophysiology.
-  Lack of caution when handling the catheter!
Injuries of the heart, perforation or tamponades.
→ Take great care when handling the catheter.
-  Improperly attached neutral electrode(s), incorrect electrodes or failure of electronic cables!
Incorrect delivery, impedance or equipment malfunction despite correct device settings
→ Ensure that the neutral electrode and its connectors are correctly attached. Avoid skin contact between different parts of the patient's body. The neutral electrode must be placed as far away from the surface ECG electrodes as possible.
→ Check the device for faults or incorrect attachment before you increase the energy supply because the power delivery is too low.

Use of unsuitable neutral electrodes!

Patient's risk of burns in the area of the neutral electrode, insufficient effect of the emitted RF power

- The neutral electrode should have the following minimum surfaces, disposable neutral electrode $\geq 136 \text{ cm}^2$, reusable neutral electrode $\geq 187 \text{ cm}^2$
- The neutral electrodes used should meet or exceed the requirements of EN 60601-2-2.

Improper handling of the catheter, generator or accessories!

Injury to the patient and/or medical personnel and/or technical facilities.

- Before using the catheter, the RF generator or other accessories, please read and understand the respective instructions for use. Follow the instructions for use.
- The catheter may only be connected to the external EP measurement point after the system has been activated.
- The connections of the cables or patient cables must be dry and clean.
- All static electricity must be kept away from the catheter system. Before connecting the catheter, the equipotential bonding cables of all devices connected to the catheter system must be connected.
- You must also take particular care to ensure that the operating table and the electronic devices which are used (e.g. X-rays) are earthed sufficiently and centrally.
- Please also note that as connecting several devices may result in a hazard due to the summation of residual currents, you must take appropriate measures to avoid this.

Conductive paths created for ablation alternating currents via the electrodes and probes of monitoring or stimulation instruments!

Risk of burns.

- To reduce the risk of burns, position the electrodes and probes as far away as possible from the ablation site and neutral electrode.

Ignition of flammable gases and other substances!

Fire or explosion.

- Ensure that highly flammable substances are removed from the working area.

The generator does not display a temperature!

Temperature control/monitoring not possible for ablation.

- Check the cable connection at the generator. If the temperature is not displayed, there may be a malfunction in the temperature sensor system. Rectify this fault before applying RF current.

Electrical contacts of the catheter or connected cables are electrically connected to the heart!

Ventricular fibrillation due to uncontrolled electrical energy

- Do not touch the plug contacts of the catheter or of connected cables with your fingers or with objects.
- Only connect the catheter's plug contacts to medical devices using connection cables which meet requirements and are intended for this purpose.
- Only use devices of protection class CF and observe their operating instructions. For defibrillation, the devices must be equipped with appropriate defibrillation protection.
- Do not allow the catheter plug or the plugs of a connecting cable to come into contact with liquids or electrically conductive surfaces. This can impair the electrical safety and functionality.

Electrostatic/electromagnetic influence on the catheter

Ventricular fibrillation due to uncontrolled electrical energy

- Only use the catheter in a properly equipped and operated medical electrophysiology laboratory (EP laboratory).

Misplacement of the catheter into the coronary system during transaortic access!

Bleeding, shock or heart attack

- Ensure that there is good imaging during transaortic access.

Long fluoroscopy times with X-rays!

Somatic and genetic damage due to high exposure to ionising radiation

- Particularly for pregnant women, children, young patients and patients who have to undergo several examinations, reconsider the risk/benefit ratio as a consequence of the potential effects of radiation.

Sterile packaging damaged or opened!

Infection of the patient due to non-sterile products

- Check the packaging before opening it and, if necessary, do not use the catheter.

Damage due to continued use contrary to the intended use and due to unsuitable storage conditions!

Damage to the central circulatory system, infection, septicaemia

- Take note of the intended use and the relevant identification on the packaging.

Improper use of the catheter or incorrect operation of the thumb switch!

Injuries, perforations and tamponades in the heart and vascular system

- The catheter may only be used as intended by an experienced cardiologist specialised in electrophysiology.
- Ensure that only trained personnel uses the catheter. Always push the thumb switch backwards to insert or remove the catheter in order to put the catheter tip into a straight position.
- If there is resistance, do not use force to push the catheter forwards or backwards.
- The catheter must not be subjected to any mechanical stress caused by manual deformation.

Reprocessing, re-sterilisation, reuse and/or technical modification of the catheter by the user or operator!

Injury to blood vessels and intracardiac injuries, infection, sepsis, death of the patient

- Vanguard Ablation Catheters uni LT G BTA/ LT BTB/ G BTC/ BTD are disposable products and exclusively intended for single use on a patient.
- The catheter may not be reprocessed, re-sterilised and reused by the user or operator.
- The catheter may not be technically modified.

Use of RF current!

Interference or destruction of implanted pacemakers and automatic implantable cardioverters/defibrillators (AICDs)

- An external pacemaker and a defibrillator must be available and ready for use throughout the entire procedure.
- Temporarily reset the pacemaker system to minimum delivery or to OFF to minimise the risk of undesired stimulation.
- Take the utmost care when the ablation point is near the atrial or ventricular electrodes. In this case, set the AICD to the OFF mode during ablation and perform a full test of the implanted equipment after

the procedure has been completed.

Ablation at a septal accessory pathway!

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- Only perform ablations at septal accessory conduction pathways after a stringent risk/benefit assessment, as the risk of a total AV block requires the implantation of a permanent pacemaker.
- An external pacemaker and a defibrillator must be available and ready for use throughout the entire procedure.

Using the catheter during Magnetic Resonance Imaging (MRI)!

Thermal-mechanical hazard possibly resulting in life-threatening internal injury to the patient

- The Vanguard Ablation Catheters uni LT G BTA/ LT BTB/ G BTC/ BTB may not be used on patients during an MRI.

Devices emitting leakage currents are operated in the immediate vicinity of the patient!

Leakage currents on catheters, lethal arrhythmias

- Ensure that no devices with leakage currents are located and operated in the immediate vicinity of the patient.

HANDLING INSTRUCTIONS

Handling the catheter requires great care to avoid injuries, perforations and tamponades to the heart and vascular system.

Storage

→Information about storage conditions is provided on the label of the product packaging.

Safety measures before use

- Before using the catheter or external pacemaker, you must read and understand the respective instructions for use.
- When connecting external devices and using them for their specific purpose, consult the corresponding instructions for use. Only use authorised medical devices under consideration of the respective instructions for use. Check all system components according to the corresponding instructions for use. The Vanguard Ablation Catheters LT G BTA/ LT BTB/ G BTC/ BTB may only be used after a successful system check.
- All static electricity must be kept away from the catheter system. You must also take particular care to ensure that the operating table and the electronic devices which are used (e.g. X-rays) are earthed sufficiently and centrally.
- It is a prerequisite for intracardiac potential recording and/or stimulation and ablation that specialist medical personnel are trained sufficiently and that a fully equipped electrophysiology laboratory is present. The procedure itself may only be carried out by a trained and experienced cardiologist specialised in electrophysiology.

Removal from packaging

- Check the sterilisation packaging for integrity and then check the sterilisation indicator. The product is sterile if the EO indicator is set to yellow. Products with damaged packaging or non-yellow EO indicator must be disposed of and a new product must be used for the procedure.
- Take the catheter out under aseptic conditions and transfer it to a sterile working environment.
- Ensure that the packaging material is disposed of completely.
- Check the catheter for obvious defects such as kinks in the shaft and loose electrodes.
- Carry out a functional check of the catheter's control mechanism and a careful visual and functional check of the connector.
- Do not use damaged catheters.

Application

- ✓ The procedure must always be performed under X-ray control and permanent ECG monitoring.
- ✓ An emergency power supply must be available to ensure continuous treatment.
- ✓ Defibrillation, emergency pacemaker therapy and acute cardioversion must be possible during the entire treatment.
- ✓ An external pacemaker, defibrillator, intubation kit, oxygen connection and emergency medication must be ready and available throughout the procedure.
- ✓ The principles of asepsis must be observed.
- ✓ In order to prevent delays in the process, a sterile substitute should be available during treatment.
- Under consideration of the aseptic rules, access a large central blood vessel through which you can insert the catheter.
- Use the appropriate connection cable to connect the catheter to a suitable recorder and/or RF generator.
- Using an X-ray and electrocardiogram as guidance, push the catheter forwards and position it in the required place in the endocardium.

NOTE: In order to make positioning the catheter tip easier, it can be bent using the thumb switch. To do so, push the thumb switch in distal direction. If you push the thumb switch in proximal direction, the tip will straighten itself again. The thumb switch may not be lubricated. The catheter is straight when the thumb switch is pushed fully into the proximal position and can then be inserted into or retracted from the port.

NOTE: The catheter may only be removed under X-ray control.

NOTE: If resistance occurs, the catheter must not be pulled back or moved with force. In this case, an X-ray must be used to check the position of the catheter.

NOTE: The port must be straightened and not moved while the catheter is being inserted or retracted. While the electrode area of the catheter is located in the insertion port, any used controllable ports may not be moved in the opposite direction of the catheter. Ports with a haemostasis valve that obstructs the catheter tip may not be used.

NOTE: In order to perform ablation under RF current, a secure connection of the catheter and the RF generator is necessary. For this purpose, only use connection cables that are compatible with the original products of VascoMed GmbH and also read the technical manual of the RF generator.

NOTE: To perform ablations, the catheter must be connected to the power output port for ablation on the RF generator. You can find more information in the Technical Manual of the relevant RF generators.

NOTE: To close the electric circuit, the neutral electrode must be connected to the reference electrode input of the RF generator.

NOTE: The multiple delivery of an RF current can lead to coagulation residues on the catheter tip. This is noticeable by the increase in impedance and temperature, and can also lead to the generator's automatic deactivation. In this case, the catheter must be removed in order to clean the electrode. In this regard, use a 0.9% physiological saline solution. Carefully clean just the tip and do not rub or bend it to avoid damaging or loosening the tip. Please also refer to the chapter on safety precautions.

NOTE: The catheters may only be used with approved recorders and RF generators that are compatible with appropriate catheters made by VascoMed GmbH. Observe the information and Technical Manual of the equipment manufacturer. Use suitable cables designed for use with appropriate catheters from VascoMed GmbH.

NOTE: For ablation, it is recommended to start with a low power and to increase it gradually as needed. A power of 75 Watt must not be exceeded during ablation. The power rating must be restricted to a maximum of 75 Watt for high power generators. When adjusting the power, take local anatomy differences in and morphology into account.

NOTE: The catheter may not be moved during ablation. Ablation in coronary vessels is prohibited.

NOTE: The signals and actual values on the display of the RF generator must be checked for coherence during the procedure. The impedance values during RF emission, the power and temperature of the head electrode can indicate catheter defects.

NOTE: The ablation temperature must not exceed 60 °C in order to ensure safe ablation with this product. The temperature displayed on the RF device is the temperature of the head electrode and not the tissue temperature.

Safety measures after use

Following the clinical procedure, the catheter must be checked for intactness for patient safety purposes. Pay attention to breaks and fragmentation of the shaft and tip. All actively implanted devices of the patient must be checked for functionality.

Disposal

Any used products must be considered as contaminated hospital waste. The used catheter is to be disposed of in accordance with the applicable national laws and directives.

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