Vanguard Diagnostic Catheter Cable
Vanguard Diagnostic Catheter Cable MTB
Vanguard Diagnostic Catheter Cable OSA

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'EMPLOI
GEBRUIKSAANWIJZING
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUÇÕES DE UTILIZAÇÃO
BRUKSANVISNING
BRUKSANVISNING
INSTRUKCJA OBSŁUGI

Deutsch
English
Français
Nederlands
Italiano
Español
Português
Norsk
Svenska
Polski

Vanguard AG
Landsberger Straße 266
12623 Berlin, Deutschland
www.vanguard.de
service@vanguard.de

C €0482

©Vanguard AG 1000002194 Revision_B_10.03.2021

Verpackungsbeschriftung und –symbole/Packaging and labelling symbols/Étiquetage et symboles sur l'emballage/ Verpakkingstekst en –symbolen/Stampigliature e simboli sull'imballo/ Etiquetas y símbolos del envase/Rotulagem e símbolos de embalagens/Påskrift og symboler på forpakningen/Förpackningsmärkning och symboler/Napisy i symbole na opakowaniu

Inhalt/Content/Contenu/Inhoud/Contenido/Conteúdo/
Contenido/Innehåll/Innhold/Zawartość/Vsebina:

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Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at service@vanguard.de or contact our field staff before using the product on patients.

General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

- ✓ Requirements
- Instructions
- Lists
-  Safety instructions

The original version of these instructions for use is in German.

Device Information

This product may only be operated if its safe use is guaranteed. Observe the safety and handling instruction!

Vanguard AG assumes no liability for damages that arise from improper operation of the Vanguard Diagnostic Catheter Cable MTB | OSA and/or disregard of these instructions for use.

Vanguard AG and the manufacturers of the original products, Medtronic, Inc. (710 Medtronic Parkway, Minneapolis, MN 55432, USA) and Osypka AG (Earl-H.-Wood-Straße 1, 79618 Rheinfelden, Germany) are not affiliated companies; the devices are remanufactured solely by Vanguard AG, who markets the device without collaboration by Medtronic, Inc. and/or Osypka AG.

“Redel” is a protected brand of the manufacturer INTERLEMO HOLDINGS S.A., Ecublens, CH or of a manufacturer associated with this company.

PRODUCT DESCRIPTION

The Vanguard Diagnostic Catheter Cable MTB | OSA are flexible, sterile connection cables intended for single use. Plug connectors, which can be designed in several standardised variants, are located at both ends of the product.

INTENDED PURPOSE

Indication

The Vanguard Diagnostic Catheter Cable MTB | OSA are intended to connect a compatible EP diagnostic catheter to a compatible EP recording system/ -stimulation device as part of an electrophysiological examination.

Guidelines for use

It is the responsibility of the physician to use a suitable surgical procedure and technology. The procedure described in the instructions for use is for information purposes only. Each physician must apply, supplement or adapt the information in the instructions for use according to his/her medical training and clinical

experience. The medical device may only be used by trained personnel in a fully equipped electrophysiology laboratory.

Contraindicationen

There are no known contraindications for the Vanguard Diagnostic Catheter Cable MTB | OSA.

TECHNICAL DATA

REF.	Corresponding REF. No. (of the original manufacturer)	Length	Connection Catheter	Connection System
35587	2ACHC (Medtronic, Inc.)	1.96 m	Redel 14-pole	10 individually protected pins
35588	990066 (Medtronic, Inc.)	1.96 m	Redel 14-pole	8 individually protected pins
35589	53820 (Osypka AG)	2.0 m	Redel 10-pole	

Other data → see label.

All Vanguard Diagnostic Catheter Cable MTB | OSA are as compatible with EP recording systems/ -stimulation devices and diagnostic catheters as the corresponding products of the original manufacturers (see table).

This product bears a CE mark in accordance with the requirements in the EC directive 93/42/EEC and all applicable changes.

SAFETY INSTRUCTIONS

The safety instructions consist of three aspects: Cause, Consequence, and Countermeasure. A safety instruction indicates a potential risk to the health or life of persons.

-  Plug and/or cable comes into contact with liquid!
Impairment with the function of the medical device
 - The plug and the cable may not come into contact with liquids.
 - The product may not come into contact with organic solvents.
-  Do not touch the plug contacts with your fingers or thin objects!
Impairment with the function of the medical device
 - Ensure that the plug contacts do not come into contact with thin objects. Do not touch the plug contacts.
-  Lack of caution when handling the Vanguard Diagnostic Catheter Cable!
Injuries to the patient, the user and/or the OR personnel

- Take great care when handling the Vanguard Diagnostic Catheter Cable.
- Observe the instructions for use of the Vanguard Diagnostic Catheter Cable.

Ignition of flammable gases or other substances!

Fire or explosion

- Ensure that flammable and easily inflammable materials are removed from the working area.

Damaged product or product is damaged during use!

Impaired electrical function of the product, misdiagnosis, injury to the patient

- Carry out a careful visual and functional inspection of the product.
- Never use damaged products.

Sterile packaging damaged or opened!

Infection of the patient due to non-sterile products

- Check the packaging before opening it. The cable may not be used if its sterile packaging is open or damaged.

Damage due to unsuitable storage conditions, use after expiry date!

Impaired electrical function of the product, misdiagnosis, injury to the patient

- Take note of the relevant identification on the packaging.
- Do not use the product once its expiry date has expired.

Reprocessing, resterilization, reuse and/or technical modification of the cable by the user or operator!

Injury, infection, death of the patient

- Vanguard Diagnostic Catheter Cable MTB | OSA are disposable products and only intended for single use on one patient.
- The Vanguard Diagnostic Catheter Cable MTB | OSA may not be reprocessed, resterilized, reused or technically modified by the user or operator.

HANDLING INSTRUCTIONS

Handling the cable requires great care to avoid injuries to the patient.

Storage

- Information about storage conditions is given on the label on the product packaging.

Safety measures before use

- The instructions for use of the Vanguard Diagnostic Catheter Cable MTB | OSA must have been read and understood before use.
- When connecting external devices and using them for specific applications, observe the corresponding instructions for use. Only use authorized medical devices under consideration of the respective instructions for use.
- Any static electricity must be kept away from the vicinity of the cable and its associated configuration.

Take particular care to ensure that the operating table and the electronic devices used (e.g. X-rays) are sufficiently and centrally earthed. Electrostatic voltages between patient and user must be equalized before the cable is used.

- A prerequisite for the use of the Vanguard Diagnostic Catheter Cable MTB | OSA is the appropriate training of medical personnel and the existence of a fully equipped electrophysiology laboratory.

Removal from packaging

- Check the sterilization packaging for integrity and then check the sterilization indicator. The product is sterile if the EO indicator is set to yellow. Products with damaged packaging or non-yellow EO indicator must be disposed of and a new product must be used for the procedure.
- Take the cable out under aseptic conditions and transfer it to a sterile working environment.
- Ensure that the packaging material is disposed of safely.
- Carry out a careful visual and functional inspection of the product. Check for defects such as kinks or damage to the plug. The product may not be used in the event of damage or malfunction.

Application

- ✓ The principles of asepsis must be observed.
 - ✓ An emergency power supply must be available to ensure continuous treatment.
 - ✓ Defibrillation, emergency pacemaker therapy and acute cardioversion must be possible during the entire treatment.
 - ✓ An external pacemaker, defibrillator, intubation kit, oxygen connection and emergency medication must be ready and available throughout the procedure.
 - ✓ The procedure must always be performed under X-ray control and permanent ECG monitoring.
 - ✓ In order to prevent delays in the process, a sterile substitute should be available during treatment.
- First connect the cable to the EP recording system/ -stimulation device and then to the diagnostic catheter, observing the marking. Check that the plug is seated correctly.
 - After completing the procedure, disconnect the cable from the catheter and the system.

NOTE: The Vanguard Diagnostic Catheter Cable MTB | OSA must be used under sterile conditions.

NOTE: If the product is damaged during use, e.g. kinks in the cable or damage to the pins in the plug, it must be replaced with a new one.

Disposal

The product may be contaminated after use and must therefore be considered as contaminated hospital waste. The cable must be disposed of in accordance with the applicable national laws and directives.

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