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GENERAL

Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at service@vanguard.de or contact our field staff before using the product on patients.

General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

- ✓ Requirements
- Instructions
- Lists
- ⚠ Safety instructions

The original version of these instructions for use is in German.

Product specifications

This product may only be operated if its safe use is guaranteed. Observe the notes on safety and usage!

VANGUARD AG assumes no liability for damages that arise from improper operation of the Vanguard Fixed Loop Diagnostic Catheter BWL and/or disregard of these instructions for use.

VANGUARD AG and the manufacturer of the original product, Biosense Webster, Inc., USA, (33 Technology Drive, Irvine, California 92618, USA) are not affiliated companies; the device is fully refurbished solely by Vanguard AG, who markets the device without collaboration by Biosense Webster, Inc.

"Biosense Webster" or "LASSO" and "LASSO 2515" are protected brands of the manufacturer of the original product Biosense Webster Inc., USA or one of its affiliated companies.

PRODUCT DESCRIPTION

The Vanguard Fixed Loop Diagnostic Catheter BWL is a circular mapping catheter with a bendable tip which has a flexible ring (loop) on its end. The loop contains several electrodes for the circular mapping of the heart's atrial structures. The catheter is brought into position in either the left or right atrium with the help of a suitable 8 F introducer sheath.

The curvature of the tip is steered by a movable plunger with a thumb switch in the hand piece at the proximal end of the catheter. The thumb switch is pushed forward to bend the tip of the shaft. When the thumb switch is in the backwards position, the tip is straight. A suitable extension cable with an appropriate

©Vanguard AG 1000001789_REV_B_09.07.2020

connector is used to connect to a standard recording device.

Intended use

The Vanguard Fixed Loop Diagnostic Catheter BWL is designed for temporary intracardiac ECG conductance and stimulation (mapping) of the atrial structures within the scope of electrophysiological examinations.

It is the responsibility of the physician to use a suitable surgical procedure and technology. The procedure described in the instructions for use is for information purposes only. Each physician must apply, supplement or adapt the information in the instructions for use according to his/her medical training and clinical experience. The medical device may only be used by a trained and experienced specialist physician.

All Vanguard Fixed Loop Diagnostic Catheter BWL have the same compatibility with accessories as the corresponding products of the original manufacturer Biosense Webster, Inc., USA.

Technical data

Shaft diameter: 7 F = 2.3 mm
Shaft diameter loop: 3 F = 1.0 mm
Usable length: 115 cm
Connection: 10 Pin (Redel)

Article No.	Original Ref No. ²	Curve type	Ring electrodes	
			Quantity	Distance [mm]
33688	D7-L10-12-RT 35O-16R	D / Lasso 12mm	10	4
33684	D7-L10-15-RT 35O-26R	D / Lasso 15mm	10	4.5
33685	D7-L10-20-RT 35O-36R	D / Lasso 20mm	10	6
33686	D7-L10-25-RT 35O-46R	D / Lasso 25mm	10	8
33690	D7-L20-15-RT 35T-26R	D / Lasso 15mm	20	4,5 (Pairs)
33689	D7-L20-20-RT 35T-36R	D / Lasso 20mm	20	6 (Pairs)
33687	D7-L20-25-RT 35T-46R	D / Lasso 25mm	20	8 (Pairs)

Other data → see label

"Redel" is a protected brand of INTERLEMO HOLDINGS S.A., Ecublens, CH or the brand of a manufacturer associated with this company.

² Corresponding Ref No. of original manufacturer
©Vanguard AG 1000001789_REV_B_09.07.2020

The catheters are used in connection with compatible connection cables from the company Biosense Webster, Inc. and suitable devices for electrophysiological mapping (recording of signals and stimulation).

This product bears a CE mark in accordance with the requirements in the EC directive 93/42/EEC and all applicable amendments.

SAFETY INSTRUCTIONS

The safety instructions consist of three aspects: cause, consequence, countermeasure.
A safety instruction indicates a potential risk to the health or life of persons.

 Proximal plug connection or handle comes into contact with liquid!

Interference with the function of the catheter

→ Prevent the connector from coming into contact with liquids.

 Electrical contacts of the catheter or connected cables are electrically connected to the heart!

Ventricular fibrillation due to uncontrolled electrical energy

- Do not touch the plug contacts of the catheter or of connected cables with your fingers or with objects.
- Only connect the catheter's plug contacts to medical devices using connection cables which meet requirements and are intended for this purpose.
- Do not allow the catheter's plug or the connection cable's plug to come into contact with liquids. This can impair the electrical safety and functionality.

 Electromagnetic interference due to an external pacemaker!

Damage or interference of other devices, ventricular fibrillation and burns

→ Only use the catheter in a properly equipped and operated medical electrophysiology laboratory.

 Plug connector damaged!

Interference with electrical safety and function

→ Carry out a careful visual inspection and functional inspection of the plug.

 Transfer of electrostatic discharges to the catheter system!

Ventricular fibrillation due to uncontrolled electrical energy

→ Only use the catheter in a properly equipped and operated medical electrophysiology laboratory.

 Long fluoroscopy times with X-rays!

Somatic and genetic damage due to high exposure to ionising radiation

→ Particularly for pregnant women, children, young patients and patients who have to undergo several examinations, reconsider the risk/benefit ratio as a consequence of the potential effects of radiation.

 Sterilised package is damaged or opened!

Infection of the patient due to unsterile products

→ Check the packaging before opening it and, if necessary, do not use the catheter.

 Damage due to improper handling and/or unsuitable storage conditions!

Damage to the central circulatory system, infection, sepsis

- Take note of the intended use and the relevant identification on the packaging.
- Do not expose the catheter to any organic solvents such as alcohol.

 Improper use of the catheter or thumb switch was not moved to the rear or loop diameter operating element was not fully turned counter-clockwise!

Injuries, perforations and tamponades in the heart and vascular system

- Ensure that the catheter is only used by trained users.
- Always push the thumb switch backwards to insert or remove the catheter in order to put the catheter tip into a straight position.
- If there is resistance, do not use force to push the catheter forwards or backwards.

 Reprocessing and reuse of the catheter by the user or operator!

Injury of blood vessels and intracardiac injuries, infection, sepsis

- Ensure that the catheter is not reprocessed.

 Improper use, turning or rotating the shaft counter-clockwise!

Entanglement of cardiac tissue in the electrode portion of the catheter

- Turn or rotate the inserted catheter only in a clockwise direction.

 Improper use, one or more ring electrodes of the Lasso-catheter in contact with the electrodes of an ablation catheter during the delivery of power!

Carbonisation on the ring electrodes

- Please ensure that the ring electrodes of the Lasso-catheter are not in contact with the ablation catheter's electrodes during energy delivery.

HANDLING INSTRUCTIONS

The catheter must be handled with great care in order to avoid injuries, perforations and tamponades in the heart and vascular system.

Storage

- Information about storage conditions is given on the label on the product packaging.

Safety measures before use

- Before using the catheter or external pacemaker, you must read and understand the respective instructions for use.
- When connecting external devices and using them for their specific purpose, consult the corresponding instructions for use. Only use authorised medical devices under consideration of the respective instructions for use.

- All static electricity must be kept away from the catheter system. You must also take particular care to ensure that the operating table and the electronic devices which are used (e.g. X-rays) are earthed sufficiently and centrally.
- A prerequisite for making an intracardiac potential recording and/or stimulation is the appropriate training of medical staff and the availability of a fully equipped electrophysiology laboratory. The procedure itself may only be carried out by a trained, experienced physician.
- The use of this catheter in pregnant woman must be considered carefully.

Removal from packaging

- Check that the packaging is not damaged.
- Take the catheter out under aseptic conditions and transfer it to a sterile working environment.
- Make sure that the packaging material is disposed of safely.
- Check the catheter for any obvious defects such as kinks in the shaft, loose electrodes or damaged plugs and, if necessary, swap it for another catheter. The catheter must not be used if there is any obvious damage.
- Carry out a careful visual inspection and functional inspection of the plug.

Application

- ✓ In order to avoid injuring the ventricular system and the heart, the intervention must always be carried out under X-ray control.
- ✓ An external pacemaker and a defibrillator must be available and ready for use throughout the entire procedure.
- ✓ The principles of asepsis must be observed.
- Apply the standard method for vessel punctures, introduction of a guide wire and the application of an introducer sheath in compliance with the instructions for use.
- Connect the catheter to a suitable recording device with an appropriate connection cable.
- Before inserting or removing the catheter, please ensure that the thumb switch used to control the bending of the tip is completely retracted. Perform and monitor the catheter movements within the patient with X-ray imaging and electrocardiograms.

NOTE: The use of the catheter is recommended together with a matching Braided Guiding Introducer Sheath PREFACE from Biosense Webster. Do not use any introducer sheaths that have holes with a diameter of more than 1.25 mm on the sides.

NOTE: Only catheters with a straight tip (retracted thumb switch) should be inserted into and removed from the introducer sheath. When inserting, ensure that the loop region is not wrinkled.

NOTE: Rotation of the catheter should only be in a clockwise direction in order to avoid entanglement with heart structure.

NOTE: In order to make positioning the catheter tip easier, it can be bent using the thumb switch. To do this, push the thumb switch forward. If you push the thumb switch backwards, the tip will straighten itself again. The catheter is straight when the thumb switch is pushed back fully and can then be removed.

NOTE: If resistance occurs, the catheter must not be pulled back with force. In this case, an X-ray must be used to check the position of the catheter.

Disposal

The catheter is to be disposed of in accordance with the applicable national laws and directives.

Contraindications

The application of the Vanguard Fixed Loop Diagnostic Catheter BWL is contraindicated for the following:

- patients with a parietal thrombus
- patients who have undergone a ventriculostomy or atriotomy in the preceding four weeks
- patients with heart valve prosthesis
- patients with active systemic infection
- patients who have shown sensitivity to foreign objects in the past or who have an allergy
- patients with severe anomalies in the heart and/or blood vessels (transposition) and obstruction

The use of transseptal procedures in patients with a thrombus or myxoma of the left atrium, interatrial variations or interatrial baffle or patch is contraindicated. Due to the risk of entanglement of the catheter in the left atrium or in the valve system, retrograde processes are contraindicated. The catheter should not be used in ventricles. Ablations may not be performed with the catheter.

Side effects

When using the catheter, breakage of the catheter, dislocation of the catheter, contact problems between the catheter and pacemaker and an increase in the stimulus threshold can cause the following undesired side effects:

- pain and discomfort in the chest
- increase in phosphokinase levels
- subcutaneous haematoma
- stimulation of skeletal muscles and nerves
- vasovagal reactions
- lesion of the brachial plexus
- lesion of the thoracic duct
- lesion of the intracardiac conduction system
- damage to heart valves and/or blood vessels (e.g. venous supplies)
- local and systemic infection
- pneumothorax due to subclavian puncture
- haemothorax
- pericarditis
- cardiac tamponade
- arteriovenous fistula
- pseudoaneurysms
- arrhythmias (e.g. ventricular tachyarrhythmias)
- cardiac thromboembolism, air embolisms
- perforation of the atrium or ventricle (myocardium/endocardium) with or without loss of simulation or ventricular fibrillation
- coronary artery spasm, coronary artery thrombosis, coronary artery transection
- transient ischaemic attack (TIA) and apoplectic stroke
- pulmonary embolism
- myocardial infarction
- total heart block
- death

Previously unknown side effects may occur at any time. The distributing company must be notified of any side effects that are not specified.

In connection with the known side effects, please also refer to the relevant specialist literature.

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