Vanguard Irrigated Ablation Catheter SJ
Vanguard Irrigated Ablation Catheter uni SJG
Vanguard Irrigated Ablation Catheter bi SJH
Vanguard Irrigated Ablation Catheter SE uni SJI
Vanguard Irrigated Ablation Catheter SE bi SJJ

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'EMPLOI
GEBRUIKSAANWIJZING
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUÇÕES DE UTILIZAÇÃO
BRUKSANVISNING
BRUKSANVISNING
INSTRUKCJA OBSŁUGI
UPUTE ZA UPOTREBU

Deutsch
English
Français
Nederlands
Italiano
Español
Português
Norsk
Svensk
Polski
Hrvatski

Vanguard AG
Landsberger Straße 266
12623 Berlin, Deutschland
www.vanguard.de
service@vanguard.de

C €2797

©Vanguard AG 1000001966_REV_C_26-08-2021

Verpackungsbeschriftung und –symbole / Packaging and labelling symbols / Étiquetage et symboles sur l'emballage / Verpakkingstekst en –symbolen / Stampigliature e simboli sull'imballo / Rotulación y símbolos del envase / Rótulo e símbolos da embalagem / Påskrift og symboler på forpakningen / Text och symboler på förpackningen / Napisy i symbole na opakowaniu / Natpisi i simboli na ambalaži

Inhalt / Content / Contenu / Inhoud / Contenuto / Contenido / Conteúdo / Innhold / Innehåll / Zawartość / Sadržaj:

Ein / one / un / een / un / un / um / én / en / jeden / jedan: Vanguard Irrigated Ablation Catheter uni SJG | bi SJH | SE uni SJI | SE bi SJJ

Medizinprodukt / Medical device / Dispositif médical / Medisch hulpmiddel / Prodotto medicale / Producto médico / Dispositivo médico / Medisinsk utstyr / Medicinteknisk produkt / Wyrób medyczny / Medicinski proizvod

X/Y
Remanufactured

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GENERAL

Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at service@vanguard.de or contact our field staff before using the product on patients.

General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

- Instructions
- Lists

Device information

This product may only be operated if its safe use is guaranteed. Observe the safety instructions, precautions and notes!


Vanguard assumes no liability for damages that arise from improper operation of the Vanguard Irrigated Ablation Catheter SJ and/or disregard of these instructions for use.


Vanguard AG and the manufacturer of the original product, St. Jude Medical (Plymouth, MN 55442 USA), are not affiliated companies; the device is fully refurbished solely by Vanguard AG, who markets the device without collaboration by St. Jude Medical. "FlexAbility", "S.JM", "Cool Point" and "Sensor Enabled" are registered trademarks of the manufacturer of the original product, St. Jude Medical (Plymouth, MN 55442 USA), or one of its affiliated companies.

"Redel" is a registered, protected brand of INTERLEMO HOLDING S.A. (INTERLEMO HOLDING AG / INTERLEMO HOLDING LTD); Chemin des Champs Courbes 28, 1024 Ecublens, Switzerland (CH).

PRODUCT DESCRIPTION

The Vanguard Irrigated Ablation Catheter SJ is intended for electrophysiological mapping and ablation in the heart. The catheter consists of a 7.5 F shaft (2) with a distal 8 F catheter tip (3) and a handle (1). The catheter tip consists of a flexible tip electrode (8) with irrigation holes (10) and three ring electrodes (9). Ablation follows with the tip electrode. This contains an integrated temperature sensor for temperature monitoring. Tip electrodes of catheters Type SJI and SJJ are additionally equipped with a magnetic sensor. This serves to visualise the position of these catheters with a compatible 3D navigation system. The electrical connection (4) for the HF generator and the 3D navigation system is proximal to the handle. The application area is flushed with saline solution through the shaft lumen via the luer connection (5) proximal to the handle. This requires the connection of a compatible irrigation pump. Potential recording and stimulation takes place with the three ring electrodes which are positioned in the working area proximal to the tip electrode. The distal catheter tip is controlled using the control elements on the handle: curvature control (6) and friction control (7). The unidirectional catheters SJG and SJI can be deflected in one direction, the bidirectional catheters SJH and SJJ in two. Correct positioning of the catheter tip is supported by the torsion-resistant shaft. The types of curvature are listed in the following table. The curvature type is marked on the label.

INTENDED PURPOSE

Indication

The Vanguard Irrigated Ablation Catheter SJ is intended for the treatment of arrhythmias in combination with compatible accessories: temporary intracardiac ECG mapping and heart stimulation, as well as for the generation of focal lesions during cardiac ablation procedures. The SJI and SJJ navigation catheters (see specifications) provide positioning information when used in conjunction with compatible S.JM™ visualisation and 3D navigation systems.

Use should be limited to patients with ventricular tachycardia where epicardial ablation is appropriate.

Guidelines for use

The medical device may only be used by a trained and experienced specialist physician. All applied medical procedures and techniques are the responsibility of the attending physician. The information contained in these instructions of use must be observed.

Contraindications

The Vanguard Irrigated Ablation Catheter SJ is contraindicated in:

- patients with active systemic infection.
- patients with prosthetic heart valves.
- patients with intracardiac thrombus or myxoma or interatrial baffle or patch via transseptal access.
- patients intolerant to anticoagulation therapy (heparin or equivalent alternative).

Side effects

A number of severe side effects have been documented for catheter ablation procedures. These include pulmonary embolism, myocardial infarction, strokes, pericardial tamponade and death.

The following complications have also been reported in previous studies or are listed in literature.

In connection with catheterisation / catheterisation procedures:

Vascular bleeding/local haematoma, thrombosis, arteriovenous (AV) fistulae, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, pericardial effusion/tamponade, thrombi, air embolism, arrhythmia and valve injuries, pneumothorax and haemothorax, water retention and lung oedema, hypoxia, pleural effusion, Acute Respiratory Distress Syndrome (ARDS), congestive cardiac insufficiency, aspiration pneumonia, pneumonia, asthma attacks, hypotension, malfunction of the implantable Cardioverter Defibrillator (ICD), anaemia, thrombocytopaenia, disseminated intravascular coagulation, haemorrhagic areas, systemic infections, urinary tract infections, apnoea caused by sedation, CO₂ retention

caused by sedation with lethargy and cholecystitis.

In connection with HF:

Chest pain/complaints, ventricular tachyarrhythmia, transitory ischaemic attack (TIA), stroke (CVA), complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery section, thromboembolism of the heart, pericarditis, perforation of the heart/pericardial tamponade, valve damage and elevated phosphokinase levels.

Independent of devices or procedures:

Urine retention, transient numbness of the limbs, Parkinson's disease and gastrointestinal haemorrhages.

Previously unknown side effects may occur at any time. The distributing company must be notified of any side effects that are not specified.

In the event of known side effects, please also refer to appropriate expert literature.

Technical data

Diameter	7.5 F (2.5 mm)
Max. diameter	8 F (2.67 mm)
Usable length	115 cm ± 1.15 cm
Irrigation port	Luer ISO 594-2 female Variant A
Sensors in catheter tip	Temperature sensor (Thermocouple) Magnetic sensor only for Type SE
Irrigation holes	lateral and distal at the tip electrode
Length of tip electrode	4 mm
Number of ring electrodes	3
Distance between electrodes	(1 – 4 – 1) mm
Rated voltage	220 V _p

SJG	uni-directional	curvature	D	---	F	---	Y
	non SE	REF No.	35096	---	35097	---	35098
		Corresponding REF no. of the OEM	A701157	---	A701158	---	A701159
SJH	bi-directional	curvature	D-D	D-F	F-F	F-J	J-J
	non SE	REF No.	35099	35100	35101	35102	35103
		Corresponding REF no. of the OEM	A701124	A701125	A701127	A701128	A701129
SJI	uni-directional	curvature	D	---	F	---	Y
	SE	REF No.	35104	---	35105	---	35106
		Corresponding REF no. of the OEM	A-FASE-D	---	A-FASE-F	---	A-FASE-J
SJJ	bi-directional	curvature	D-D	D-F	F-F	F-J	J-J

SE	REF No.	35107	35108	35109	35110	35111
	Corresponding REF no. of the OEM	A-FASE-DD	A-FASE-DF	A-FASE-FF	A-FASE-FJ	A-FASE-JJ

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All Vanguard Irrigated Ablation Catheter SJ are as compatible with the accessories (e.g. generator, cable) as the corresponding products of the original manufacturer St. Jude Medical (see table). When the catheter is connected to an HF generator that provides an HF output voltage $> 220 V_p$, the maximum HF output voltage must be set to $220 V_p$.

SAFETY INFORMATION

1. Cardiac ablation procedures may only be performed by physicians who have mastered the techniques of high-frequency catheter ablation. These procedures must be performed in a fully equipped electrophysiological laboratory.
2. Only approved medical devices may be used during the procedure with strict adherence to their instructions for use.
3. The temperature values transmitted by the sensor in this catheter are not the actual temperatures of the tissue. The temperature values only represent the values of the irrigated electrode.
4. Due to the duration and intensity of X-rays for fluoroscopic imaging during catheter ablation procedures, there is a risk of significant exposure to radiation with acute radiation damage and an increased risk of somatic and genetic effects for the patient and laboratory staff. A catheter ablation procedure may only be performed after thorough consideration of the potential radiation exposure and after the implementation of adequate radiation protection measures. The long-term risks of prolonged exposure to X-rays were not determined. For pregnant women and pre-pubescent children, a catheter ablation procedure must therefore be carefully considered.
5. Pacemakers and implantable cardioverters/defibrillators can be impaired or damaged by high-frequency signals. Discharges from intracardiac defibrillators (ICD) during application may result in injury to the patient. The following precautions should therefore be taken during use: a) Provide external stimulation and defibrillation during ablation; b) Deactivate intracardiac defibrillators (ICD) as they may be discharged and injure the patient or damage the patient during ablation; c) Take special care during ablation near permanent stimulation or defibrillation electrodes; and d) Carefully inspect implantable aggregates in all patients following ablation.
6. Special care must be taken when placing lesions near the native conduction system.
7. The long-term risks of HF ablation lesions have not been investigated, especially with regard to lesions near the native conduction system.
8. Ablation in and near arterial coronary vessels has been associated with myocardial infarction and death.
9. To avoid overhydration, the patient's fluid balance should be monitored throughout the application. Proceed according to your hospital protocol.
10. Before inserting the catheter into the vascular system, be sure to remove all air from the tubes and

catheter as trapped air can cause injury or death.

11. The materials of the catheter are not compatible with magnetic resonance imaging (MRI).
12. The catheter may get caught in heart or blood vessels, which is a possible complication in electrophysiological applications.
13. There is a risk of vascular perforation or dissection with each placement of the electrodes. To avoid heart damage, perforation or tamponade, move the catheter with caution.
14. Do not apply excessive force if you encounter resistance when advancing or retracting the catheter.
15. If you are using an electrophysiological measuring station, it must either be isolated from the patient or have an isolated patient cable.
16. The right phrenic nerve can be injured during ablation along the vena cava.
17. This product is intended for single use only. It may not be reprocessed or reused. Reuse may result in product failure, patient injury and/or transmission of infectious diseases from one patient to another patient or third party.
18. The safety and efficacy of epicardial ablations have not yet been adequately investigated.

PRECAUTIONS

1. The HF energy is switched off in the event of an alarm signal from the Cool Point™ irrigation pump. Communication and liquid flow must be checked. Check the irrigation line for blockages such as kinks and air bubbles. Remove the catheter from the patient and check both the catheter and the electrodes. If necessary, clean the electrodes with sterile saline solution and gauze. Before reinserting the catheter, ensure that the irrigation ports are unblocked and irrigate the catheter.
2. During HF ablation, the impedance display of the HF generator must be monitored continuously. If the impedance suddenly increases without exceeding the preset limit value or if steam suddenly escapes, HF ablation must be stopped manually. The doctor must assess the situation. If necessary, the catheter must be removed from the patient and any clots removed from the distal tip. If the catheter is defective, it must be replaced by a new catheter. Before reinserting the catheter, ensure that the irrigation ports are unblocked. After completing the test, reinsert the catheter for another HF ablation test.
3. At power levels above 40 watts, the sudden escape of steam may occur frequently. At the maximum permissible power level (50 watts), increased collateral damage may also occur. This power level should only be selected if the desired result cannot be achieved with a lower power level.
4. The catheter must be advanced under fluoroscopic imaging control. This minimises the risk of heart damage, perforation or tamponade. Compatible catheter positioning and navigation systems can be used in conjunction with fluoroscopic imaging.
5. Before inserting and withdrawing the catheter, always straighten the curvature using the control elements.
6. Always ensure continuous irrigation to avoid coagulation in and around the electrodes.
7. Excessive bending or buckling can damage the catheter. Manual pre-bending of the distal end may result in damage to the control mechanism and injury to the patient.
8. It has been proven that ablation systems with irrigation produce larger lesions than HF ablation catheters without irrigation. Special care should be taken when ablating near structures that are sensitive to electrical charges and near thin-walled or other arterial structures.

9. The use of suitable filters of the mapping systems is necessary to enable continuous monitoring of surface ECG or IEGM during the application of high-frequency energy. Monitoring systems with built-in devices for limiting the HF current are recommended for this purpose.
10. Needle monitoring electrodes are not recommended.
11. Do not immerse the handle and cable plug in liquid. This can impair the electrical function.
12. The connecting cables may not be in contact with the patient and other electrical discharges.
13. If the irrigation flow is interrupted, the catheter must be checked immediately outside the patient and irrigated again. Before reinserting the catheter into the patient, the irrigation flow must be ensured.
14. Do not ablate without using the Cool Point™ irrigation pump.
15. The tip electrode may not be twisted or pulled. Excessive force may detach the electrode from the catheter shaft.
16. After use, dispose of the product and packaging in accordance with your hospital regulations or observe official and/or local guidelines.
17. The catheter must not be brought into contact with organic solvents such as alcohol.

HANDLING INSTRUCTIONS

Storage

The device must be stored in dry and light-protected conditions at 5°C to 25°C. Information about storage conditions is given on the label on the sterile packaging.

Removal from packaging

- Check that the packaging is not damaged and check the validity of the expiry date. Do not use the device if the packaging is damaged or opened and the expiry date has expired.
- Remove the catheter from its packaging under aseptic conditions and transfer it to a sterile working environment. Do not throw the instrument.
- Check the catheter for damage such as kinks in the shaft, loose electrodes, damaged plugs or defects such as clogged irrigation ports, difficulty in adjusting the distal end and insufficient fixation of the distal end. Do not use damaged or inadequate catheters.

Application

- Ensure that the generator and the necessary accessories are set up according to the diagram in the HF generator user manual. Do not connect the compatible cable until the catheter has been connected and prepared as described below. Ensure that all unused terminals of the 4-electrode IEGM cable are isolated when used with a mapping or recording system. This reduces the risk of unwanted current paths to the heart.
- Connect a sterile saline-filled luer-lock syringe to the luer port on the catheter. Push the contents of the syringe into the catheter to ensure that all irrigation ports are unblocked.
- Connect the catheter with the (standard type) luer connections to the irrigation system. It must be possible to start up the pump at a minimum injection pressure of 1 bar.
- Ensure that all air bubbles have been removed from the tubing and catheter prior to insertion. The catheter must be irrigated with a high flow rate setting. Add heparin to the saline infusion solution in accordance with the patient's anticoagulation state.

- Connect the compatible connecting cable to the catheter while observing the connection polarity. Do not use force when connecting the plug as this could damage the contacts. Connect the cable to the connector socket labelled "Catheter Extension Connector" on the front of the generator. For information on connecting to compatible magnetic tracking systems, please refer to the cable's user manual.
- Switch on the generator and the pump. Information on the operation, installation and communication of both devices can be found in the HF generator's user manual.

Pump settings

Minimum continuous flow rate	120 ml/h (2 ml/min)
Flow rate during ablation	600 - 999 ml/h (10 - 17 ml/min)

HF Generator settings

Output	Recommended initial setting	20 W
	Maximum setting	50 W
Temperature	Recommended initial setting	40 C
	Maximum setting	45 C
Duration	Maximum ablation at one site (uninterrupted)	60 s

- Before selecting the ablation parameters on the generator, ensure that the indifferent electrode is correctly attached to the patient.
- Select 20 watts as the initial power.
- Select 40°C as the initial temperature.
- Ensure that the catheter tip is in the neutral (straightened) position before insertion. For easier insertion, you can use an introducer sheath with a minimum diameter of 8.5F. To prevent irrigation line occlusion, the catheter must be irrigated continuously while it is in the vascular system. Continue with irrigation until the catheter has been removed from the body.
- Connect the catheter to the compatible catheter positioning and navigation systems using a compatible cable. For information on the use of catheter positioning and navigation systems in conjunction with this device, refer to the relevant instructions for use.
- Under consideration of the aseptic rules, access a large central blood vessel through which you can insert the catheter.
 - The catheter should be inserted through a peripheral vessel and its position monitored under fluoroscopic imaging.
 - The distal catheter tip is deflected using the control elements on the handle. In the unidirectional catheter, the curvature is set by feeding and retracting the curvature control element. In the bidirectional catheter, the curvature is set by turning the curvature control element. Do not pull on the luer connection.
- Observe the impedance display of the HF generator before, during and after HF ablation. If the impedance suddenly increases during HF ablation without exceeding the preset limit, HF ablation must

be stopped manually. The doctor must assess the situation. If necessary, the catheter must be removed from the patient and any clots removed from the distal tip.

- To start the HF therapy (ablation), press the "RF Energy Delivery" button on the generator. The irrigation flow rate is automatically increased by the pump from the initial continuous flow rate to a high flow rate.
- For setting a drag lesion, move the catheter linearly, whereby you may not remain on one site for longer than 60 seconds.
- Starting with the recommended initial power of 20 watts, the power can be increased to a maximum of 50 watts as needed to achieve an effective lesion. IEGM and impedance must be assessed before changing the power setting. At a power level above 40 watts, the sudden escape of steam may occur more frequently. The power should only be increased to these values if the desired result cannot be achieved with lower power levels.
- If the output temperature of 40 °C is reached without the preset output power having been reached, the temperature may be increased to the permissible maximum (45°C). Check IEGM and impedance beforehand.
- In the event of suddenly escaping steam, the HF energy output must be interrupted.
 - After the sudden escape of steam or automatic deactivation, remove the catheter and visually check for clots, charring or other catheter defects.
- If the pump issues an alarm and the irrigation process is interrupted, the catheter must be immediately removed from the patient, inspected and irrigated again (see generator and pump user manuals). At the end of each ablation, the pump automatically returns to the initial continuous flow rate with the selected delay.
- At the end of the procedure, the catheter must be placed in the neutral (straightened) position before it can be removed from the patient.

NOTE: In the bidirectional catheter, the degree of curvature fixation can be adjusted all the way to a full lock. This is set with the friction control element. The degree of fixation increases the further the knob is turned clockwise (in plus direction). Delivery follows in non-fixed setting.

NOTE: In the unidirectional catheter, the resistance of the control element for curvature adjustment and therefore also the degree of fixation of the curvature can be adjusted up to the locking point. In the delivery state, the elements are set to the highest resistance or the highest degree of fixation. If an adjustment is required, the friction control element can be rotated. To reduce the resistance or the degree of fixation, rotate to the right and to the left to increase it.

Disposal

The product must be disposed of in accordance with the applicable national laws and directives.

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