## VANGUARD Reference Patches UNI

- EN** INSTRUCTIONS FOR USE  
VANGUARD Reference Patches UNI suitable for CARTO® 3
- DE** GEBRAUCHSANWEISUNG  
VANGUARD Reference Patches UNI geeignet für CARTO® 3
- FR** MODE D'EMPLOI  
VANGUARD Patches de référence UNI convient pour CARTO® 3
- IT** ISTRUZIONI PER L'USO  
VANGUARD Elettrodo di riferimento UNI adatte per CARTO® 3
- NL** GEBRUIKSAANWIJZING  
VANGUARD Referentiepatches UNI geschikt voor CARTO® 3
- NO** BRUKSANVISNING  
VANGUARD Referanseplastre UNI egnet for CARTO® 3
- PL** INSTRUKCJA OBSŁUGI  
VANGUARD Elektrody odniesienia UNI przeznaczone do systemu CARTO® 3
- PT** INSTRUÇÕES DE UTILIZAÇÃO  
VANGUARD Patches de Referência UNI adequados para CARTO® 3
- ES** MANUAL DE INSTRUCCIONES  
VANGUARD Parches de referencia UNI aptos para CARTO® 3
- SV** BRUKSANVISNING  
VANGUARD Referenspatchar UNI lämpliga för CARTO® 3

VANGUARD AG  
Landsberger Straße 222  
12623 Berlin, Deutschland  
[www.vanguard.de](http://www.vanguard.de)  
[service@vanguard.de](mailto:service@vanguard.de)

## Symbol Definitions / Symboldefinitionen / Définition des symboles / Definizione dei simboli / Symboldefinities / Symboldefinisjoner / Definições dos símbolos / Definiciones de los símbolos / Symboldefinitioner

EN

Content / Inhalt / Contenu / Contenuto / Inhoud / Innhold / Zawartość / Conteúdo / Contenido / Innehåll

DE

Reference Patches / Referenzelektroden / Patches de référence / Elettrodo di riferimento / Referentiepatches / Referanseplastre / Elektrody odniesienia / Patches de Referência / Parches de referencia / Referenspatchar

FR

Instructions for use / Gebrauchsanweisung / Mode d'emploi / Istruzioni per l'uso / Gebruiksaanwijzing / Bruksanvisning / Instrukcja obsługi / Instruções de utilização / Manual de instrucciones / Bruksanvisning

IT

For single use / Zur einmaligen Verwendung / Prévu pour une utilisation unique / Monouso / Bestemd voor eenmalig gebruik / Til engangsbruk / Do jednorazowego użytku / Para utilização única / De un solo uso / För engångsbruk

NL

Do **NOT** use if packaging is damaged / Bei beschädigter Verpackung **NICHT** verwenden / NE PAS UTILISER si l'emballage est endommagé / Se la confezione non è integra **NON** utilizzare / **NIET** gebruiken als de verpakking beschadigd is. / Må **IKKE** brukes dersom forpakningen er skadet / **NIE** używać, jeżeli opakowanie jest uszkodzone / **NÃO** utilizar se a embalagem estiver danificada / **NO** utilizar en caso de que el envase esté dañado / Produkten får **INTE** användas om förpackningen är skadad

NO

PL

Expiry date / Verwendbar bis / À utiliser avant le / Utilizzare entro / Uiterste gebruiksdatum / Siste bruksdato / Nie używać po tej dacie / Data limite de utilização / Utilizable hasta / Sista förbrukningsdag

PT

Item number / Artikelnummer / Numéro d'article / Codice articolo / Artikelnummer / Artikkelnummer / Nr artykułu / Referência / N.º de artículo / Artikelnummer

ES

Lot number / Losnummer / Numéro de lot / Numero di lotto / Batchcode / Batch-kode / Kod partii / Código de lote / Código de lote / Partikod

SV

Serial number / Serien-Nummer / Numéro de série / Numero di serie / Serienummer / Serienummer / Numer seryjny / Número de série / Número de serie / Serienummer

Store away from light / Lichtgeschützt aufbewahren / Conserver à l'abri de la lumière / Conservare la riparo dalla luce / Beschermd tegen licht bewaren / Oppbevares beskyttet mot lys / Chronić przed światłem / Conservar ao abrigo da luz / Almacenar protegido de la luz / Förvaras ljusskyddat

Store in a dry place / Trocken aufbewahren / Conserver au sec / Conservare in luogo asciutto / Droog bewaren / Oppbevares tørt / Przechowywać w suchym miejscu / Conservar em local seco / Consérvese en un lugar seco / Förvaras torrt

Information about acceptable ambient temperature during storage and transportation / Angaben zur zulässigen Umgebungstemperatur bei der Lagerung und beim Transport / Indications relatives à la température ambiante autorisée pour le stockage et le transport / Indicazioni relative alla temperatura ambiente ammessa per l'immagazzinamento e il trasporto / Specificaties voor de toelaatbare omgevingstemperatuur bij opslag en transport / Angivelser om tillatt omgivelsestemperatur ved oppbevaring og transport / Dane dotyczące dopuszczalnej temperatury otoczenia w czasie przechowywania i transportu / Informações sobre a temperatura ambiente admissível durante o armazenamento e o transporte / Indicación de la temperatura ambiente para el almacenamiento y el transporte / Uppgifter om tillåten omgivningstemperatur vid förvaring och transport

Follow the instructions for use / Gebrauchsanweisung beachten / Respecter le mode d'emploi / Consultare le istruzioni per l'uso / Houd u aan de gebruiksaanwijzing / Følg bruksanvisningen / Przestrzegać instrukcji obsługi / Observar as Instruções de Utilização / Seguir el manual de instrucciones / Beakta bruksanvisningen

Electronic instructions for use / Elektronische Gebrauchsanweisung / Mode d'emploi électronique / Istruzioni d'uso elettroniche / Elektronische gebruiksaanwijzing / Elektronisk bruksanvisning / Elektroniczna instrukcja obsługi / Instruções de uso eletrônicas / Instrucciones electrónicas de uso / Elektronisk bruksanvisning

Non-sterile / Unsteril / Non stérile / Non sterile / Niet steriel / Ikke steril / Niesterylne / Não esterilizado / No esterilizado / Ikke-steril

Manufacturer / Hersteller / Fabricant / Produttore / Fabrikant / Produsent / Producent / Fabricante / Fabricante / Tillverkare

CE mark / CE Zeichen / Label CE / Marchio CE / CE-markering / CE-merke / Znak CE / Marcação CE / Marca CE / CE-märkning

Medical device / Medizinprodukt / Dispositif medical / Prodotto medicale / Medisch hulpmiddel / Medisinsk utstyr / Wyrób medyczny / Dispositivo medico / Producto medico / Medicinteknisk produkt

Manufacturing date / Herstellungsdatum / Date de fabrication / Data di produzione / Fabricagedatum / Produksjonsdato / Data produkcji / Data de fabric / Fecha de fabricación / Tillverkningsdatum

EN

DE

FR

IT

NL

NO


PL

PT


ES

SV


1


2

3

4

## ENGLISH

### GENERAL

Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at [service@vanguard.de](mailto:service@vanguard.de) or contact our field staff before using the product on patients.

#### General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

-  Safety instructions
- Instructions
- Lists

The original version of these instructions for use is in German.

#### Device information

This product may only be operated if its safe use is guaranteed. Observe the notes and safety precautions! VANGUARD AG assumes no liability for damages that arise from improper operation of the VANGUARD Reference Patches UNI and/or disregard of these instructions for use. "CARTO® 3" is a protected brand of the medical devices manufacturer Biosense Webster, Inc., 3333 Diamond Canyon Road, Diamond Bar, CA 91765, USA.

### PRODUCT DESCRIPTION

The "VANGUARD Reference Patches UNI" are surface electrodes with holders for magnetic sensors for use in combination with a "CARTO® 3 System" from Biosense Webster, Inc. USA. All Vanguard Reference Patches are as compatible with accessories as the corresponding products from the original manufacturer Biosense Webster, Inc. The electrodes are self-adhesive. A plastic bracket for magnetic sensors is attached in the centre on the rear. The electrode has an 80 cm long connection cable. A connection socket for the sensor unit is located at the end of the cable.

#### Intended purpose

The "Vanguard Reference Patches UNI" are intended for

connection to the "CARTO® 3 System" from Biosense Webster, Inc. for non-fluoroscopic catheter navigation. They serve to attach the magnetic sensors to the patient and, together with the "CARTO® 3 System", enable the transmission of the catheter's position.

It is the responsibility of the physician to use an appropriate surgical procedure and technique. The procedure described in the instructions for use is for information purposes only. Each physician must apply, supplement or adapt the information in the instructions for use according to his/her medical training and clinical experience. The "VANGUARD Reference Patches UNI" may only be used by appropriately trained and experienced cardiologists and trained clinical surgical staff specialised in electrophysiology.

#### Declaration of conformity

 This product bears a CE mark in accordance with the requirements in the REGULATION (EU) 2017/745.

### SAFETY INSTRUCTIONS

The safety instructions consist of three aspects: Cause, Consequence, Countermeasure. A safety instruction indicates a potential risk to the health or life of persons.

-  Application on damaged skin!  
**Occurrence of skin intolerances on the patient**
  - Do not use the patches on damaged skin. In the event of skin intolerances, terminate the application immediately.
-  Application on insufficiently prepared skin!  
**The patches detach, which causes faulty signal transmission, which results in extended surgery and/or in incorrect measurement and/or incorrect diagnosis/treatment**
  - Only use the patches on sufficiently prepared skin.
-  Contact with organic solvents such as alcohol!  
**Detaching electrodes, impaired patch functionality**
  - Avoid contact with organic solvents. Leave the organic solvents applied to the patient to dry before positioning the patches.
-  Positioning the patches over pacemakers, ICDs or comparable devices!  
**Risk of burns, impaired patch or device functionality**
  - Patches may not be applied over pacemakers, ICDs or comparable devices.

Item No.	Corresponding ref. no. of the original manufacturer	Electrode diameter	Cable length	Internal socket diameter
35802	CREFP6	8 cm	80 cm	2 mm

- ⚠ Contact or overlapping with ECG patches, neutral electrodes or other reference sensors!  
**Risk of burns, impaired patch functionality**  
→ Avoid contact with other patches or electrodes.
- ⚠ Using the patches on other devices than those intended!  
**Risk of burns, impaired patch or device functionality**  
→ VANGUARD Reference Patches may only be used with the CARTO® 3 system.
- ⚠ Positioning the patches over bony areas!  
**Interference with electrical safety and function**  
→ Do not position the patches over bony areas such as the scapula or sternum.
- ⚠ Tensile stress on the cable!  
**Interference with electrical safety and function**  
**Risk of damages to the sensor unit**  
→ Carefully remove the electrodes from the carrier film. Do not pull on the cable.  
→ When separating the cable from the sensor unit, hold the plug connection and loosen the connection with a gentle pulling and rotating movement.  
→ Avoid pulling on the cable.
- ⚠ Pulling the cable/sensor!  
**The patches detach and/or the sensors are damaged, which causes faulty signal transmission, which results in extended surgery and/or in incorrect measurement and/or incorrect diagnosis/treatment.**  
→ Avoid pulling on the sensor. Check the sensor's firm seating.  
→ Open the holders to remove the sensors. Do NOT pull on the sensor cable.  
→ Avoid pulling on the cable.
- ⚠ Cable is in contact with liquid!  
**Impaired patch function**  
→ Prevent the cable from coming into contact with liquids.
- ⚠ Cable/Patches damaged!  
**Interference with electrical safety and function**  
→ Carry out a careful inspection of the patches.
- ⚠ Plug connection is loose!  
**Interference with electrical safety and function**  
→ Check all plug connections for perfect fit.
- ⚠ Storage of unused patches!  
**Faulty patch function due to unsuitable storage**  
→ Store unused patches in their sealed, original packaging. Patches not stored in their sealed, original packaging may not be used. Information about storage conditions is provided on the label on the product packaging.
- ⚠ Multiple use!  
**Risk of cross-contamination with injury, illness or death**  
→ The product is intended for single use on one patient.

- ⚠ Using together with reference sensors from other manufacturers!  
**Incorrect diagnosis, impaired treatment success**  
→ Do not use VANGUARD Reference Patches together with reference sensors from other manufacturers.
- ⚠ Magnetic sensors insufficiently padded!  
**Risk of pressure sores and damages to the patient**  
→ Place the patient on a sufficiently soft support.

## HANDLING INSTRUCTIONS

- The product is not sterile.
- Only intended for single use on one patient.
- In the event of skin intolerance, terminate application immediately.

### Storage

- Information about storage conditions is given on the label on the product packaging.
- Do not use the patches after their expiry date (📅).
- Store unused patches in their sealed, original packaging.

### Safety measures before use

- Before using the patches, you must read and understand the respective instructions for use.
- When connecting external devices and using them for their specific purpose, consult the corresponding instructions for use. Only use authorised medical devices under consideration of the respective instructions for use.
- All static electricity must be kept away. You must also take particular care to ensure that the operating table and the electronic devices which are used (e.g. X-rays) are earthed sufficiently and centrally.

### Removal from packaging

- Check that the packaging is not damaged. If the product packaging is damaged, dispose of the product and use a new product for surgery.
- Ensure that the packaging material is disposed of safely.
- Check the patches for obvious faults such as damaged cables and replace the patches if necessary. The patches may not be used if there are any obvious damages.

### Application

- After removing the patches from their packaging, check them for any damage.
- Position the patches around the heart and inside the CARTO® 3 mapping zone.
- Positioning on the back:
  - connect (3) patches with the sensor unit and attach the sensors in the appropriate brackets, refer to Figure 1.
  - Then position the patches on the back of the seated patient as follows:

- localise the centre of the heart slightly to the left of the spine at the same level as the seventh vertebra,
  - attach the first patch in the centre approximately 12 cm above the centre of the heart,
  - attach the second patch 12 cm to the right below the centre of the heart,
  - attach the third patch 12 cm to the left below the centre of the heart.
  - The patches on the back should form a triangle - refer to Figure 3.
  - Once the patches are attached to the back, bring the patient into a reclining position.
- Place the location pad (positioning base) under the patient table so that the patches on the patient's back are in the centre of the location pad. Use "Location Setup" in the CARTO® 3 system to check the position of the location pad. If required, correct the position of the location pad.
- Positioning on the chest:
- connect (3) patches with the sensor unit and attach the sensors in the appropriate brackets, refer to Figure 1.
  - Then position the patches on the chest of the reclining patient as follows:
    - attach the first patch on the right lower ribcage (5th/6th intercostal space) on the right central collarbone line.  
Use "Location Setup" to check precise positioning inside the mapping zone.
    - attach the first patch on the left lower ribcage (5th/6th intercostal space) on the frontal axillary line (below the V5 ECG patch).  
Use "Location Setup" to check precise positioning inside the mapping zone.
    - Attach the third patch underneath the collarbone, to the left or right of the sternum (depending on the position of the defibrillator electrodes or whether an ICD or pacemaker must be taken into consideration). The patches on the chest should form a triangle - refer to Figure 4. Use "Location Setup" to check precise positioning inside the mapping zone.
- Use "Location Setup" in the CARTO® 3 system from Biosense Webster to check that all chest patches are located inside the mapping zone.
- After use, remove the sensors from the brackets and disconnect the patches from the sensor unit.
- To remove the sensors from the patches, open the brackets manually, see Figure 2. Do not pull on the sensor cable.
  - To separate the cable from the sensor unit, hold the plug connection and loosen the connection with a gentle pulling and rotating movement.

NOTE: The patches are intended for use in combination with the "CARTO® 3 System" from Biosense Webster, Inc.

NOTE: Before use, ensure that the patient's skin surface is completely shaved and dry.

NOTE: Also pay attention to the safety information.

NOTE: If a patch cannot be correctly attached to the patient's skin, discard and use a new one.

Any used products must be considered as contaminated hospital waste. The used product must be disposed of in accordance with the applicable national laws and directives.

## Contraindications

There are no known contraindications regarding the use of the "VANGUARD Reference Patches UNI".

## Side effects

The following undesirable side effects may occur when using "VANGUARD Reference Patches UNI":

- Reddened skin.

Previously unknown side effects may occur at any time. The distributing company must be notified of any side effects that are not specified.

In connection with the known side effects, please also refer to the relevant specialist literature.

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ANVISNING: Elektroderna är avsedda att användas i kombination med CARTO® 3 systemet från Biosense Webster, Inc.